Hartmetallinstrumente

Tungsten Carbide Instruments | Instrumentos de Carburo Tungsteno

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MEISINGER tungsten carbide instruments impress with their particularly stable and functional design, their precise, fast cutting toothing for optimized work, outstanding performance with a long service life and optimal concentricity. The instruments are made of high-density fine-grain carbide.

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63

TWIST FINIERERTWIST FINISHER
TWIST FINISHER

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

Crown cutters are one of the most important instruments in the dental practice today. The variety of different restorative materials offers an equally wide variety of different requirements for the crown cutting. A quick, easy and, not least, economic separation of crowns and bridges requires high-quality instruments that are specially adapted to the properties of the restorative materials in order to make removal particularly efficient. Therefore, MEISINGER offers the optimum crown cutter for every material. All instruments are specialists in their respective field of application and convince with high stability, best cutting performance and longevity.

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Instructions for use:

- Cutting angle of about 45°
- Optimum speed of about 160.000 rpm
- Use with the red contra-angle with low contact pressure

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Crown Cutters with Cross Cut

Materials:

- Hard precious metal and non-precious metal alloys

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






Crown Cutters with Diamond Cut

Materials:

- Especially hard metals, non-precious metal alloys and veneering ceramic

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Diamond Crown Cutters




Materials:

- Extreme hard monolithic ceramic, zirconia

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




- Greatly reduces cut vibration
- Optimal revolving properties
- Fast excavation as a result of efficient cutting geometry
- Minimal feeding pressure
- Easy recognition due to blue color ring

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
The MEISINGER Black Cobra tungsten carbide burs are characterized by their outstanding performance characteristics. Due to their extreme hardness, the instruments have a very high resistance to breakage and reduced wear. The sharp precision cutting edges with aggressive cross-cut toothing ensure the lowest friction values and, as a result, particularly efficient and quiet instruments.

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
The following surface examination underlines the special usability of the MEISINGER Twist Finishing Burs in comparison to standard Tungsten Carbide and Diamond instruments. For this, the instruments have been applied on composite. Then, the roughness of the treated surface has been measured.

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The examination shows that the MEISINGER Twist Finishing Burs are highly suitable for processing of composite materials due to their precise twist angle design. The surface and achieve perfect, extremely smooth results. The cutting and non-grinding material processing also achieves a smoother surface than ultrafine diamond instruments.

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- The different cutting edges ensure a perfect working result and an extremely smooth surface
- Very smooth running
- Highly suitable for composites and other materials
- Different color rings clearly differentiate the instruments

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Note:
Particularly well suited for removing adhesive residues in the bracket zone!

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


Plain cut, medium (Ref.-No. 175), blue ring, for effective trimming of hard materials and for obtaining especially smooth surfaces

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
Plain cut, coarse (Ref.-No. 215), green ring, for rapid and effective trimming of acrylics and plaster

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
X-cut, extra fine (Ref.-No. 141), yellow ring, for fine, precise trimming of ceramics, precious metals and all non-precious alloys

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

X-cut, fine (Ref.-No. 140), red ring, for delicate and difficult trimming of hard alloys and acrylics

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
“MINIS” like FX, x-cut, fine (Ref.-No. 140), red ring, for fine and accurate trimming of hard alloys and acrylics

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X-cut, medium (Ref.-No. 190), blue ring, for trimming large surfaces of non-precious alloys and acrylics

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
“MINIS” like GX, x-cut, medium (Ref.-No. 220), blue ring, for trimming non-precious alloys and acrylics

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
X-cut, coarse (Ref.-No. 220), green ring, for trimming large surfaces of acrylics, tray materials, and plaster

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
X-cut, super coarse (Ref.-No. 224), black ring, the optimum for acrylic, plaster, and wet plaster

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
Special toothing for non-precious metal alloys, coarse (Ref.-No. 137), green ring

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
Special toothing for PEEK, fine (Ref.-No.146), red/grey Ring

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
Special toothing for the processing of thermoplastic materials, blue Ring, for the finishing of rails

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
Cross cut, fine (Ref.-No. 134), red ring, for fine trimming on all alloys and hard acrylics

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
Cross cut, medium (Ref.-No. 176), blue ring, for effective and smooth trimming on all prosthesis and surgical acrylics

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
Special X-cut (Ref.-No. 191), blue ring, for use on model casting alloys, acrylics, and plaster

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
- Reinforced, resistant, and widened beveled cutters
- Ample chip spaces for the expulsion of material
- Specially suited for all permanent acrylics, acrylic resins, and gypsum material
- In spite of high material removal, a good surface quality is reached
- High rotation speeds, very high output (up to 30,000 min⁻¹)

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



The excellent features of the Black Cobra carbide cutters are enhanced by the special CARBOCER® (diamond-like carbon) coating. A specially developed process applies a diamond-like carbon coating to the Cobra carbide cutters with extremely hard, wear-resistant, low-friction, corrosion-resistant surfaces while retaining the original surface properties. This surface refinement leads to protection against aggressive cleaning agents. Contamination is also easier to detect and remove.

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Carbide chip coated for effective trimming on soft reline material. Reduced heat generation during processing due to the large surface area of the working parts. The cooling effect of the Z4 is additionally enhanced by the hollow shape.

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
MEISINGER Flexible Polishers guarantee odorless, precise, vibration-free performance in any case, from challenging ultra-fine grinding to high gloss polishing. Their use guarantees optimal surface quality on all materials. The optimal adjustment to the different material groups ensures the lowest heat development during processing. The instruments consist of high-quality fine abrasives and polishing agents in a resilient, flexible and latex-free compound of polysiloxane and synthetic rubber mounted on stainless steel shaft.

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
2-Step Systems, Single-Use Diamond Composite Polishers

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


MEISINGER Twist Polishers were specially developed for efficient high gloss polishing. The flexible polishing lamellae allow adaptation to any surface structure and produce flawless results even in areas that are difficult to access, such as the occlusal surface. Due to its material-friendly design, the original surface structure of the restoration is retained. The different polishing stages also ensure short processing times and extremely long service lives.

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Application:

- The Twist Polishers should always be used with very light pressure.
- To protect the polishing lamellae, use over sharp margins should always be avoided.
- Please note the recommended speeds.

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Intraoral application of Twist Polishers for the polishing of composites and hybrid ceramics:

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Prophy Twist Polishers

Prophy Twist Polishers / Prophy Twist Polishers


Empfohlene Drehzahl in min ⁻¹ : Recommended Speed in rpm: Revoluciones recomendadas en min ⁻¹ :	Maximale Drehzahl in min ⁻¹ : Maximum Speed in rpm: Revoluciones máximas en min ⁻¹ :
7.000 - 12.000	20.000

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Universal Polishing system for:

Gold alloys
Amalgam
Acrylics
Cement

- 1. Step
Pre-polishing
- 2. Step
Polishing
- 3. Step
High shine polishing

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Universal Polishing system for:

Precious metals
Non-precious metals
Model casting
Amalgam

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Universal Polishing system (unmounted) for:
Precious metals
Acrylics

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
PROPHYLAXEPOLIERER | POLISHERS FOR PROPHYLAXIS | PULIDORES PROFILÁCTICOS

Prophy Twist Polishers | Prophy Twist Polishers | Prophy Twist Polishers

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- Disposable instruments with plastic shank
- Flexible polishers lamellae allow adaptation to any surface structure
- Especially for polishing regions with difficult access
- Perfect results within the shortest time

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- Disposable instruments with plastic shank
- Polishing into gingival pockets
- Easy and cost-efficient use
- Application without polishing paste

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Prophy-Strips | Prophy-Strips | Prophy-Strips

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Blau / Rot
Grob / Medium

Gelb / Weiß
Fein / Extra Fein

The strips are made of high flexible polyester and are one-side coated with aluminium oxide particles.

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Prophy-Clean | Prophy-Clean | Prophy-Clean

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Particularly robust polishers made of white corundum in a plastic bond. **Scale Removal, precise removal of tooth discolorations, easy cleaning of areas difficult to reach, no damages to enamel**

- coarse grain
- medium grain

- Removal of tartar, plaque, and discolorations
- Removal of residues of bracket bonding
- Removal of excess cement
- Subgingival root smoothing

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Polishers for prophylaxis, Cup
(With polishing paste)

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The wide collar of the polisher protects the contra-angle handpiece optimally against the penetration of the polishing paste.

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Flexible polishing discs were designed **for working with composites, ceramics and metals.**

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Achieving the "high shine" on your restorations

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Especially for simple polishing of composite restorations

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Note: To be used with little contact pressure! |

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Diamond-polisher for:
soft composites

- 1. Step**
Pre-polishing
- 2. Step**
Polishing
- 3. Step**
High shine polishing

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Polishing system:**For polishing of composites**

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Diamond polishing system:

Especially for hard, highly-filled composites of the latest generation and nano composites.

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
Twist Polisher:

Especially for simple polishing of all composite restorations, also in difficult to access regions.

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Occlu Twist Polisher:

- Especially for occlusal polishing of composite restorations
- Reaches even the deepest fissures
- Gentle on material
- High-gloss results

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Diamond polishers for:**hybrid composites**

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**Diamond impregnated porcelain
polishing system:**
Ceramics

- 1. Step**
Pre-polishing (super coarse)
- 2. Step**
Pre-polishing (coarse)
- 3. Step**
Polishing (medium)
- 4. Step**
High shine polishing (fine)
- 5. Step**
Mirror finish polishing (extra fine)

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
Diamond impregnated porcelain polishing system:

Especially for simple polishing of all ceramic restorations (except hybrid ceramics), also in regions difficult to access.

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2. Step
Polishing (medium)

3. Step
High shine polishing (fine)

4. Step
Mirror finish polishing (extra fine)

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Diamond impregnated polishing system for:

Silicate ceramics
Leucit-reinforced ceramics
Lithium Silicate (e.max®)

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Diamond impregnated polishing system for:
CAD/CAM restorations
Oxide ceramics
Zirconium oxide

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**Porcelain polishers for
adjusting, smoothing, and polishing of:
Ceramics**

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Model casting polishers for:
Non precious metals
Model casting

Brown
Very hard and abrasion resistant

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Polishers for acrylics for:

Acrylic Teeth, PMMA, Hybrid Ceramic, PEEK

- Particularly suitable for the polishing of provisional supplies
- Complete polishing process in just one step
- Elevated pressure for surface preparation
- Reduced pressure for high shine polishing

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PMMA1 - PMMA/PEEK Adjusting
and Polishing Kit

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Silicone polishers with aluminum oxide for:

Synthetic materials for prosthetics

Polishers are distinguished by a long service life

- 1. Step
Pre-polishing (coarse)
- 2. Step
Polishing (fine)

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Polishing system for: Synthetic materials for prosthetics

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Cleaning instructions:

- Soap suds, no aggressive disinfectants
- Ultrasonic application possible
- No resistance to thermal disinfection or sterilisation

Recommended Speed: 10.000 rpm
Maximum Speed: 15.000 rpm

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The MEISINGER LUSTER® Paste is a 1-step diamond polishing paste used for high-shine polishing of all dental ceramic materials. It is suitable both for extraoral and intraoral use and gives a perfect high-shine to all ceramic restorations. The high-quality polishing paste significantly simplifies the polishing of occlusal surfaces and the maintenance of existing restorations. The MEISINGER LUSTER® Paste contains diamond particles and is available with two different grain sizes: 8-12 µm (fine) and 16-18 µm (medium).

Intended for use on all ceramic materials

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The MEISINGER LUSTER® Zirconia Paste is a 1-step diamond polishing paste used for extraoral high-shine polishing of zirconia restorations. The high-quality polishing paste significantly simplifies the polishing of occlusal surfaces. The MEISINGER LUSTER® Zirconia Paste contains diamond particles with a grain size of 4-8 µm.

Intended for extraoral use on zirconia

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MEISINGER ceramic abrasives for practice and laboratory are made of high-quality corundum and silicon carbide. They impress with their uniform grain size and optimal bonding hardness and are characterized by a high-strength, stainless steel shank. At the same time, they are optimized in shape and size for high-precision, intraoral material abrasion. The high production quality guarantees a precise and central concentricity.

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Easy processing with low heat development

The MEISINGER diamond porcelain grinders are suitable for easy processing of zirconium oxide, Lithium disilicate, feldspathic ceramics and veneering ceramics. The two different grains, coarse and medium, allow good removal with very little heat development. They are very well suited for thinning the crown margins. Additionally, the diamond porcelain grinders can be used without water cooling.

Especially for the treatment of zirconium oxide restorations.

- synthetic, diamond-coated binding
- good removal
- low heat development
- no water cooling necessary
- good reduction of crown margins

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
Diamond grinding tools made of silicon carbide for the safe machining of sensitive zirconium oxide surfaces.

- dry machining
- no microfissures
- Application without pressure

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Diamond abrasives for precise and safe processing of ceramics.

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White „Arkansas“ abrasives, extra fine, for smooth grinding and polishing of filling materials and composites.
Significant reduction of the edge gap.

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Corundum abrasives (Al_2O_3), pink, for use on all types of metals.

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Abrasives “K+B”, blue – aluminium oxide for NON Precious Alloy

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Mounted abrasives for chrome cobalt alloys, white, hard bonding

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- Low speed range (recommended: 5.000-7.000 rpm)
- Minimized heat development due to opened-pored structure
- Quick and efficient
- Multifunctional applications
- Very long-lasting
- No melting

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Use the following unmounted discs with Mandrel 305RF, Size 50
Information on page 146

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MEISINGER steel instruments are manufactured from one piece tungsten-vanadium steel or stainless steel. They demonstrate quality through their functional design, precision toothing, and above average concentricity. The instruments offer high elasticity and optimal material hardness.

HSS-Steel instruments are produced from high-strength steel, and therefore, offer an additional heat strength as well as high endurance.

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Martensitic, stainless steel, shaped for surgical applications

- ⇒ High corrosion resistance is achieved with a higher amount of chromium than with tungsten vanadium steel
- ⇒ Long service life for surgical or similar applications

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Selected tungsten-vanadium steel quality

- ⇒ Low corrosion resistance
- ⇒ Long service life with regard to hard materials
- Sterilization just possible with suitable pretreatment as well as testing for cleaning and disinfection unit suitability.
- Approx. 15% greater hardness than stainless steel instruments

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Selected high performance steels (High Speed Steel)

- ⇒ Not corrosion resistant, but higher share of chromium than with tungsten vanadium instruments
- ⇒ Long service life also with regard to very hard-materials

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For oral and maxilla-facial surgery MEISINGER offers one of the most extensive lines of allport burs, surgical cutters and trephines. The optimized construction of the instruments guarantees safe and fast achievement of desired results.

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

After tooth extraction or a longer period of edentulism, unwanted unevenness, rough surfaces or sharp edges may develop on the jawbone. In this case, the MEISINGER Ridge Contouring Bur offers the possibility of surgical contouring and modification of the alveolar ridge to smooth the bone. So, the alveolar ridge is prepared optimally for implant insertion or denture restoration.

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The special easy-cutting and efficient toothing of these instruments, as well as the particular twist angles, allows for a wide range of application of the instruments within the field of oral surgery.

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- Conservative surgical bone treatment (i.e. apicectomies, exposure and cutting of impacted teeth, preparation of a bone flap as well as general osteotomy in bone substance)
- Cutting of hard tooth substance, in particular tooth roots
- Axial penetration of bone is possible (z.B. in apicectomies)
- Easy-cutting toothing
- Narrow neck offers better view of the field of work
- Working parts with reduced diameter for narrow osteotomies
- Extremely low vibration due to optimized smooth running properties

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
¹ Only for RF. Pre-treatment is necessary for instruments made of tool steel. See the notes on preparation.

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









Special care should be exercised when using trephines. Specifically, the recommended speeds are not to be exceeded.
Recommended speeds: 300 - 500 rpm

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



! Technology by Dr. Bernd Giesenhausen

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The MEISINGER gingiva trimmers are intended for use in dental mucosal surgery. They have been specially developed for the versatile and gentle treatment of gingival tissue. The working part consists of particularly high-quality and stable zirconia. The main application focus of the MEISINGER gingiva trimmers is the safe and gentle removal of unwanted granulation tissue and hyperplastic gingiva (papillectomy) – even in hard-to-reach areas. Furthermore, the trimmers are used in the enlargement of the sulcus after crown preparation, for the best possible exposition of the preparation margin and the creation of an exact impression. The gingiva trimmers are also particularly suitable for the tissue-saving exposure of covered healed implants and for use in periodontal treatment. The application takes place without cooling, so that the resulting rotational energy can be used for tissue modulation. Thanks to heat coagulation, bleeding tendency is reduced.

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

Instructions for use:

- Application without cooling
- Optimal speed: 300.000 - 450.000 rpm

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For Ear, Nose, and Throat Surgery, MEISINGER offers a multitude of different special instruments made of various materials (steel, tungsten carbide, diamond). To guarantee the specific application in surgery, the instruments are available in a shank length from 70 mm to 125 mm.

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HNO-Bohrer | E.N.T. burs | Fresas O.R.L.

INNENGEKÜHLTES INSTRUMENT | INTERNALLY COOLED INSTRUMENT |

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Stainless steel surgical cutters with internal cooling system according to Kirschner. Special handpieces necessary (e.g. NSK, KaVo, W&H, Micro-Mega).

Note: Products equipped with a conduit to supply a cooling medium and other products with lumens (canals, drillings etc.):

- Actively rinse the lumen during pre-cleaning, manual cleaning and the disinfection process


- Products containing lumens without feed-through channels must not be re-used

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MEISINGER special instruments and product lines for endodontics offer the dentist a wide range of instruments for the safe and simple implementation of all common endodontic treatment methods.

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COMPREHENSIVE ACCESS KIT	
COMPREHENSIVE ACCESS KIT	
COMPREHENSIVE ACCESS KIT	169

MEITRAC I ENDO SICHERHEITSSYSTEM	
MEITRAC I ENDO SAFETY-SYSTEM	
MEITRAC I PREVENCIÓN DE ENDODONCIA	169

WURZELKANALINSTRUMENTE | ROOT CANAL INSTRUMENTS | ENSANCHADORES DE CONDUCTOS

176


Fig.	Schaft / Shank	Ref.-Nr. / Ref.-No.		6	6	6	6	6	6
			L mm	26,0	26,0	26,0	26,0	26,0	26,0
				0,15	0,20	0,25	0,30	0,35	0,40
176		340 814 655 453		015	020	025	030	035	040

Feilen nach Hedstroem, langer Griff | Hedstroem files with long handle | Limas según Hedstroem con mango largo

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Sortiment, Art-Nr. 17610 (015-040)
Assortment, Art-No. 17610 (015-040)
Surtido, ref. 17610 (015-040)

180GR


Fig.	Schaft / Shank	Ref.-Nr. / Ref.-No.		10	10	10	10	10	10
			L mm	19,0	19,0	19,0	19,0	19,0	19,0
180GR	RA XL	330 206 679 336		050	070	090	110	130	150

Wurzelkanalweiterer "Gates", rostfrei | Root canal reamers "Gates", stainless | Ensanchadores radiculares, inoxidables "Gates"

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Sortiment, Art-Nr. 18001 (050-150)
Assortment, Art-No. 18001 (050-150)
Surtido, ref. 18001 (050-150)

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

- Special instruments for root canal preparation
- Long narrow necks allow for an excellent view onto the working field
- Optimal, low-vibration concentricity
- Highly-stable construction
- Less dust generation than ultrasonic tips
- Color code allows for quick identification and selection of the instruments

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Recommended Speed: 500 - 800 rpm
Maximum Speed: 2.000 rpm

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For an optimized process of coronal root canal preparation, the Endo Access Bur Kit provides you with a compilation of all highly stable MEISINGER Endo Access Burs:

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Art.-No. ENDO2

COMPREHENSIVE ACCESS KIT | COMPREHENSIVE ACCESS KIT | COMPREHENSIVE ACCESS KIT

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Safe access into the root canal is imperative for any successful root canal treatment. That is why MEISINGER created the Comprehensive Access Kit, a carefully curated selection of diamond and tungsten carbide instruments necessary for a safe and predictable endodontic access preparation. This kit provides all the burs needed to obtain straight line access into the canal through different materials such as zirconia, PFMs, and enamel. Conveniently organized in an autoclavable bur block, this kit is essential for any root canal treatment.

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
With this endodontic safety-system, broken fragments can be grasped and removed successfully, precisely and quickly, (e.g. root canal instruments, root filling pins, silver points, and root-posts from the root canal). This system is designed in such a way that fragments with diameters of 0.15 to 0.50 mm can be removed safely.

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MEISINGER offers a variety of instruments that support many different aspects and working needs in the area of Orthodontics.

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Prophy-Clean	
Prophy-Clean	
Prophy-Clean	173

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For removing adhesive residues and finishing the tooth surface

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Particularly robust polishers made of white corundum in a plastic bond.

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For finishing plastic, plaster and NEM alloys

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





For the processing of deep-drawn rails, cooling grooves avoid clogging of the working part

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


Whether preparation, polishing or targeted processing of CAD/CAM restorations, MEISINGER offers perfectly coordinated system solutions in pre-assembled product lines. The very time-consuming compilation of instrument sets is therefore no longer necessary for the user. Due to the close cooperation with practitioners and leading scientists, who co-develop and test the systems, they are optimally adapted to the everyday requirements.

When it comes to the optimum polishing of CAD/CAM restorations, MEISINGER is the first choice. The individually tailored kits for each common restoration material are specially adapted to the different requirements of e.g. zirconium oxide, lithium disilicate, hybrid ceramics or zirconium-reinforced lithium disilicate. For all materials MEISINGER offers both a laboratory set for extraoral work and a surgery set for intraoral work. The kits consist of single or multi-stage polishing systems, which can be used without polishing paste and which have proven to be good by having excellent high gloss and perfect adaptation to the natural tooth surface. In addition to the long-life diamond polishers, many systems contain other instruments for grinding contact points or contact surfaces. The instruments are well organised in the high-quality drill stand and can be conveniently sterilised after use.

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
PRÄPARATION PREPARATION PREPARACIÓN	180-188	LUS07 LUSTER® Extraoral Twist Kit Large	191	LUS02 LUSTER® for Hybrid Ceramics Intraoral Adjusting and Polishing Set	200	2638 Acrylic Adjustment Kit	211
DDA01 CEREC® Kit Digital Dental Academy	180	JK03 LUSTER® CAD/CAM Lab Kit	192	GC05 Twist Polisher Lab Kit for Hybrid Ceramics according to Greg Campbell, DDS	201	VAL01 Flexible Acrylic Adjusting and Polishing Kit	212
2565 Occlusal Reduction Kit developed with Dr. Athas N. Kometas, DMD, PA	181	SC01 CAD/CAM Finishing Kit	192	GC06 Intraoral Twist Polisher Kit for Hybrid Ceramics according to Greg Campbell, DDS	201	TITAN TITANIUM TITANIO	212
2565T Tapered Occlusal Reduction Kit	181	SILIKAT KERAMIKEN SILICATE CERAMICS CERÁMICAS DE SILICATO	193-194	2663 Twist Polishing Kit	202	1321 Titanium Power Kit developed with University RWTH Aachen	212
SC02 CAD/CAM Prep Kit	182	LUS80 LUSTER® for Lithium Disilicate Adjusting and Polishing Kit	193	2571 Twist Finishing Kit	203	ENDO ENDO ENDO	213-214
BC02 Black Cobra Crown Preparation Kit	182	LUS100 LUSTER® for Lithium Silicate Intraoral Adjusting Kit	193	LUS66 LUSTER® Composite Finishing Kit	203	ENDO1 Comprehensive Access Kit	213
JK02 Prep and Finishing Kit according to James Klim, DDS	183	GC01 Twist Polisher Kit for Silicate Ceramics according to Greg Campbell, DDS	194	2664 Composite Finishing Kit	204	ENDO2 Endo Access Bur Kit	213
2561 Athen Preparation Crown & Inlay Kit according to Stavros Pelekanos, DDS	183	GC02 Intraoral Twist Polisher Kit for Silicate Ceramics according to Greg Campbell, DDS	194	SF10 Super Flexible Discs (10mm)	204	2271 MEITRAC I	214
2450 Athen Preparation Veneer Kit according to Stavros Pelekanos, DDS	184	ZIRKONOXID ZIRCONIA CIRCONIO	195-198	SF14 Super Flexible Discs (14mm)	205	PERIIMPLANTITIS PERIIMPLANTITIS PERIIMPLANTITIS	214
2595 Opti-Shape Set according to Dr. Gabriele Diedrichs	184	LUS85 LUSTER® for Zirconia Adjusting and Polishing Kit	195	FELDSPATKERAMIKEN FELDSPATH CERAMICS CERÁMICAS DE FELDESPATO	206-209	2575 Periimplantitis Kit	214
2522 Preparation Set according to Prof. Dr. Karl-Peter Meschke	185	LUS91 LUSTER® for Zirconia Intraoral Adjusting Kit	195	LUS41 LUSTER® Lab Porcelain Polishing Kit	206	SONSTIGES MISCELLANEOUS OTRO	215-217
2525 Microdia Kit	185	GC03 Twist Polisher Lab Kit for Zirconia according to Greg Campbell, DDS	196	LUS30 LUSTER® for Porcelain Intraoral Polishing Kit	206	BC04 Black Cobra Complete Crown Removal Kit	215
2413 Laminate Veneer Kit	186	GC04 Intraoral Twist Polisher Kit for Zirconia according to Greg Campbell, DDS	196	JK01 LUSTER® Intraoral Polishing Kit according to James Klim, DDS	207	2740 Black Cobra Cutter Set	215
DW001 Discover White Veneer Kit by Dr. Mariana Mintcheva	186	GC07 Adjusting and Polishing for Pre & Post Sintered Zirconia Extraoral Kit according to Greg Campbell, DDS	197	2609 CAD/CAM Polishing Kit	207	2540 Plaque Set	216
2531 Preparation Set according to Dr. Andres Baltzer	187	2670 Zirkon Master Kit	197	LUS60 LUSTER® CAD/CAM Porcelain Preparation and Polishing Kit	208	2530 Fuzzies Assortment Kit	216
2685 Cercon® 0° / 2°	187	5901 Zirkon Diamond Kit	198	1470 / 1471 Ceramics Polishing Kit	208	JK05 Acryl / Resin Universal Kit developed with Dr. James Klim, DDS	217
DB01 CAD/CAM Prep Kit Anterior according to Dr. Daniel Buttermann, DDS	188	JK04 Zirconia Finishing Kit according to Dr. James Klim, DDS	198	2605 CAD/CAM Advanced Kit	209		
DB02 CAD/CAM Prep Kit Posterior according to Dr. Daniel Buttermann, DDS	188	ZIRKONVERSTÄRKTES LITHIUM SILIKAT ZIRCONIA REINFORCED LITHIUM SILICATE SILICATO DE LITIO REFORZADO CON CIRCONIO	199	2662 Intra Oral Porcelain Adjustment Kit	209		
ALL CERAMICS ALL CERAMICS ALL CERAMICS	189-191	LUS03 LUSTER® for Zirconium reinforced Lithium Silicate Extraoral HP Laboratory Set	199	KUNSTSTOFFE ACRYLICS ACRÍLICO	210-212		
CKK01 CHARD ADVANCED DENTISTRY CAD/CAM CERAMIC CHARACTERISATION KIT according to Dr. Simon Chard	189	LUS04 LUSTER® for Zirconium reinforced Lithium Silicate Intraoral Adjusting and Polishing Set	199	PMMA1 PMMA/PEEK Adjusting and Polishing Kit	210		
LUS06 LUSTER® Intraoral Twist Kit	190	HYBRIDKERAMIKEN / KOMPOSITE POLY CERAMICS / COMPOSITE CERÁMICAS HÍBRIDAS / COMPOSITE	200-205	BC03 Black Cobra Crown Temporization Kit	210		
LUS05 LUSTER® Extraoral Twist Kit	191	LUS01 LUSTER® for Hybrid Ceramics Extraoral HP Laboratory Set	200	2582 Temporary Solution Kit developed with ZTM Rainer Michel	211		

DDA01 CEREC® Kit Digital Dental Academy

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The MEISINGER CEREC® Kit has been developed in cooperation with international trainers at the Digital Dental Academy. This kit lets you experience significant simplification and time savings in your daily work. The attractive feature of the CEREC® Kit is its logical structure that guides you through your whole work procedure: Excavation – Preparation – Polish. Along with the most important instruments for the precise preparation of inlays, onlays, partial crowns and crowns, the kit also includes the outstanding MEISINGER crown cutters HM7RX for metals and Z838L for zirconia. The high-quality diamond polishers for a flawless finish complete the contents.

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The Patented Occlusal Reduction Burs (United States Patent number 6,511,322) allow for precise occlusal reductions to be achieved. The vertical orientation provides great access to the occlusal plane. The Occlusal Reduction Burs are available in 1.5 mm, 1.8 mm, 2.0 mm, 2.2 mm, and 2.4 mm depths corresponding to proper occlusal reduction for various restorative materials. A 1.0 mm bur is also available for labial reduction. Developed by Athas N. Kometas, DMD.

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The Tapered Occlusal Reduction Diamonds have an 8° taper and a unique self-limiting angled collar. The Tapered Occlusal Reduction Diamonds are excellent for producing angled walls for inlay and onlay preparations, as well as standard occlusal depth reduction. The Tapered Occlusal Reduction Diamonds are available in 6 color-coded depths (1.0mm, 1.5mm, 1.8mm, 2.0mm, 2.2mm, and 2.4mm) and have two bands to signify the difference from the Original Occlusal Reduction Diamond Burs.

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The CAD/CAM Prep Kit contains an ideal selection of award-winning diamonds including all the necessary shapes for inlays, onlays and full coverage crowns. Included in the kit are MEISINGER's 1.5 mm and 2 mm depth limiting Occlusal Reduction Burs for predictable occlusal reductions.

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BC02 Black Cobra Crown Preparation Kit

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MEISINGER's BC02 kit is a multi-purpose crown preparation kit that carries our exclusive Black Cobra diamonds for quick, efficient tooth reduction and includes matching trimming and finishing carbides to define and smooth the preparation.

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According to Dr. James Klim, DDS, Preparation architecture, smooth flowing margins and line angles are the cornerstones for successful milled indirect restorations. Dr. Klim's Prep and Finishing Kit includes Black Cobra diamonds for efficient shaping and complementary carbide finishing burs for finishing the margins. The kit also includes MEISINGER's patented occlusal reduction burs for true occlusal reduction everytime.

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2561 Athen Preparation Crown & Inlay Kit

according to Stavros Pelekanos, DDS

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The Athens Preparation & Crown Inlay Kit provides an intelligent compilation of instruments for optimal crown and inlay preparation. The torpedo-like diamond instruments equipped with an axially extending locator pin are suited for the creation of a defined chamfer. The circular chamfer is easily prepared supragingivally, while the polished locator pin is guided around the tooth without pressure.

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The Athen Preparation Veneer Kit comprises an optimal selection for veneer application in the anterior region. At the same time, the 834 groove burs facilitate the preparation of a defined cutting depth on the labial surface of the teeth. The 852 diamond instruments are particularly conducive to the rapid removal of the hard dental tissue (labially) down to the predetermined depth, whereby the instrument can also be used palatally / lingually.

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2595 Opti-Shape Set

according to Dr. Gabriele Diedrichs

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This assortment for preparation for crowns, partial crowns, inlays and onlays is especially designed for galvano-ceramic, metal-ceramic and all-ceramic restorations.

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2522 Preparation Set

according to Prof. Dr. Karl-Peter Meschke

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An easy and user friendly preparation set for defined material removal for crowns and inlay preparation.

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The minimally invasive diamond alternatives to oscillating instruments.

In the field of microdentistry, the microdia kit offers the possibility of obtaining healthy tooth substance to the maximum extent possible during preparation. The rotating diamond instruments are distinguished by small working parts and extremely slim, fitted instrument shafts. They provide a very good view of the working field and are therefore ideally suited for preparation technology. The highest precision and select diamond grit of the instruments facilitate the performance of high-definition work with a low level of contact pressure, especially in the case of microscope preparations.

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The two depth cutters 834-016 and 834-021 assist with creating pre defined depths on facial surfaces. The smaller creates a 0.3 mm depth on the central incisors and canines, while the larger head size creates a 0.5 mm depth on lateral incisors. The 852 series reduces the facial tooth structure to the pre determined depth. The 852 shaped diamond may also be used on the lingual wrap. The 820 shaped diamond removes interproximal flash and overhang. The 889L and 898 diamonds create an infinite margin and may assist subgingivally.

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The Discover White Veneer Kit, developed by Dr. med. dent. Mariana Mintcheva, contains specially adapted instruments for veneer preparation. Included is the MEISINGER depth marker with which a specific depth of preparation can be defined. Furthermore, the kit contains the innovative MEISINGER CF-Diamond, which has two different grain sizes and is especially suitable for the preparation of veneers. With the finely covered tip, a perfect preparation margin can be created without changing the instrument. The Super Flexible Discs in three different grits finally ensure the perfect finish of the restoration.

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This preparation set enables the dentist, first of all, to achieve optimal material strength with pre-defined cutting during the chamfer and step preparation for all-ceramic without special aids. The torpedo-like diamond instruments with an axially extending locator pin are suitable for preparing a defined chamfer. The circular chamfer is easily prepared supragingivally, while the polished locator pin is guided around the tooth without pressure.

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2685 Cercon® 0° / 2°

Cercon® is a registered trademark of Dentsply Sirona, Germany

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Assortment for Cercon® in CAD/CAM-Technology. Instruments used for the parallel milling technique have an angle of 0° and for the conical milling technique an angle of 2°.

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¹ 314=FG

DB01 CAD/CAM Prep Kit Anterior

according to Dr. Daniel Buttermann, DDS

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
The MEISINGER CAD/CAM Prep Kit Anterior was developed in cooperation with Dr. Daniel Buttermann, DDS, and is ideally suited for the preparation of CAD/CAM restorations in the anterior region. For an optimal design of the preparation margin, the kit contains matched coarse and fine diamond instruments. This optimises your workflow and helps you to work as efficiently and consistently as possible. The DB01 Kit contains depth markers for even reduction in the anterior region, thus providing the optimum basis for highly esthetic restorations.

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The MEISINGER CAD/CAM Prep Kit Posterior was developed in cooperation with Dr. Daniel Buttermann, DDS, and is ideally suited for the preparation of CAD/CAM restorations in the posterior region. For an optimal design of the preparation margin, the kit contains matched coarse and fine diamond instruments. This optimises your workflow and helps you to work as efficiently and consistently as possible. The DB02 Kit contains the patented MEISINGER „Occlusal Reduction Burs“, with which a defined reduction of the occlusal surface can be achieved for glass-ceramic and zirconium oxide restorations. With the selected, reduced-length instruments, ideal axial reductions can be achieved even in the posterior region.

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¹ 314=FG

ALL CERAMICS | ALL CERAMICS | ALL CERAMICS

CCK01 Chard Advanced Dentistry CAD/CAM Ceramic Characterisation Kit

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Chard Advanced Dentistry CAD/CAM Ceramic Characterisation Kit - Designed for characterisation of chairside ceramic restorations including Empress, emax® and hybrid restorations. This kit has been designed with the aesthetic digital dentist in mind to give creative control to surface texture, lustre and anatomy after the restoration is milled. This allows the dentist to provide predictably beautiful results if they are just polishing or if they are staining and glazing.

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The LUSTER® Intraoral Twist Kit has been especially developed for intraoral high-gloss polishing of all ceramic restorations. The flexible polishing discs make it possible to adapt to any surface structure and also render flawless results in areas that are difficult to reach, such as occlusal surfaces. Due to designs that are gentle on materials, the original surface structure of the restoration is maintained. In addition, the grades of diamond grit separated into four parts provide for short processing times and extremely high levels of durability.

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
The LUSTER® Extraoral Twist Kit has been especially developed for the simple high-gloss polishing of all ceramic restorations. The flexible polishing discs make it possible to adapt to any surface structure and also render flawless results in areas that are difficult to reach, such as occlusal surfaces. Due to designs that are gentle on materials, the original surface structure of the restoration is maintained. In addition, the grades of diamond grit separated into four parts provide for short processing times and extremely high levels of durability.

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

LUS07 LUSTER® Extraoral Twist Kit Large

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According to James Klim, DDS, "The combination of laboratory quality shaping diamonds and diamond polishers are designed for efficient and safe CAD/CAM ceramic finishing. This kit is what I use for finishing my CEREC® restorations in my clinical dental theater and hands-on teaching classes".

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This universal contouring and finishing kit was created to provide a wide range of solutions for today's ceramic restorations. The kit provides all of the necessary tools to produce excellent results for both anterior and posterior crowns. Included in the kit are diamonds and abrasives for sprue removal and contouring of line angles, including adjusting and contouring anatomy. Our newest Twist Polishers for Ceramics bring out the luster of your restorations entirely without paste while removing minimal material and safe enough to use over the margins without damage. The instruments are labeled and organized for easy reference in an anodized bur block.

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¹ 104=HP

SILIKATKERAMIKEN | SILICATE CERAMICS | CERÁMICAS DE SILICATO

LUS80 LUSTER® for Lithium Disilicate Adjusting and Polishing Kit

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This kit features the MEISINGER specially designed polishers for Lithium Silicate restorations. The 9736G wheel allows for adjustments, and the two-step polishing system will re-glaze the restoration.

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This kit contains specially formulated diamonds to make adjustments to Lithium Silicate as well as the two-step Lithium Silicate polishing system to re-glaze Lithium Silicate intraorally after the restoration has been placed.

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
¹ 204=RA, 314=FG

GC01 Twist Polisher Kit for Silicate Ceramics

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Developed in collaboration with Greg Campbell, DDS, the MEISINGER Twist Polisher Kit for Silicate Ceramics is designed to be simple and efficient. Included in the kit are the necessary tools to adjust and polish Silicate Ceramics including an abrasive and diamond for sprue removal and contouring anatomy and a sintered (layered) diamond S805-018 for enhancing occlusal anatomy. MEISINGER's new Twist Polishers for Ceramics bring out the luster of your restorations and safe enough to use over the margins without damage. This entire complement of tools can be used on pre-crystallized and crystallized restorations with remarkable results.

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GC02 Intraoral Twist Polisher Kit for Silicate Ceramics

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Developed in collaboration with Greg Campbell, DDS, the MEISINGER Twist Polisher Kit for Silicate Ceramics is designed to be simple and efficient. This kit provides the proper diamonds and polishers for adjusting and polishing Silicate Ceramics. The diamonds will gently adjust the occlusal surfaces and the margins in preparation for polishing. Two Polishing points and MEISINGER's Twist polishers are included for blending and polishing the margins, cusps, and occlusal surfaces.

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LUS85 LUSTER® for Zirconia

Adjusting and Polishing Kit

BruxZir® is a registered trademark of Glidewell Laboratories

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The MEISINGER LUSTER® for Zirconia has been specially developed for use on layered and monolithic Zirconium Oxide restorations. The hardness of Zirconia materials ranges from 850 MPa to 1200 MPa. One of the keys to adjusting and polishing these materials is to perform the tasks quickly and without creating excessive heat. MEISINGER has formulated these zirconia polishers to minimize heat buildup and thus minimize the possibilities of micro-fractures without excessive wear and tear of the polishing materials. The kit contains two-step polishers that can be used without polishing paste. Further instruments for adjusting and contouring are included.

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GC03 Twist Polisher Lab Kit for Zirconia

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Developed in collaboration with Greg Campbell, DDS, the Twist Polisher Kit for Zirconia is designed to be simple and efficient. Included in the kit are only the necessary tools to adjust and polish Zirconia and aid in creating beautiful and healthy smiles. The abrasives and diamonds are perfect for sprue removal and contouring anatomy. A sintered (layered) diamond S807-018 is included for long service life for enhancing occlusal anatomy. MEISINGER's new Twist Polishers for Zirconia bring out the luster of the restoration while removing minimal material. These Twist polishers are safe enough to use over the margins without damage.

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GC04 Intraoral Twist Polisher Kit for Zirconia

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Developed in collaboration with Greg Campbell, DDS, the Twist Polisher Kit for Zirconia is designed to be simple and efficient. Included in the kit are only the necessary tools to adjust and polish Zirconia and aid in creating beautiful and healthy smiles. The Intraoral kit contains the perfect diamonds and polishers for the adjustment and polishing of Zirconia. The polishing points and Twist polishers polish and blend the restoration to produce the perfect luster.

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¹ 314=FG, 204=RA

GC07 Adjusting and Polishing for Pre & Post Sintered Zirconia Extraoral Kit

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Developed in collaboration with Greg Campbell, DDS, the Adjusting and Polishing for Pre & Post Sintered Zirconia Extraoral Kit for Zirconia is designed to be fast, repeatable, and predictable. Dr. Campbell has specifically selected instruments to adjust and polish pre-sintered Zirconia as well as post-sintered Zirconia and aid in creating beautiful and healthy smiles. The kit contains the necessary tools to adjust and polish both pre & post sintered Zirconia.

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
The Zirkon Master Kit for use on zirconium oxide restorations. These high-tech materials give pleasure to an ever growing popularity and allow the new zirconium abrasives to be used easily and precisely. Minimal grinding pressure makes this an ideal instrument for designing in CAD/CAM-technology with minimal heat buildup. The special composition of the grinders significantly reduces the formation of microcracks even during dry processing.

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The instruments of the Zirkon Diamond Kit have been specifically developed for working on Zirconium Dioxide surfaces. Various grain sizes and a number of different shapes offer flexible and efficient use for different indications. The diamond grains are well adapted to adjust Zirconium Oxide Ceramics.

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The Zirconia Finishing Kit was developed by Dr. James Klim, DDS, and is a combination of tungsten carbide instruments, an abrasive, sintered diamonds and polishers. These selected instruments are ideal for pre-sintered zirconia and the finishing of zirconia restorations. All instruments are suitable for use in the laboratory and can be conveniently stored in a high quality bur block.

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The MEISINGER LUSTER® for Zirconium reinforced lithium silicate polishing system was specially developed for use on zirconium reinforced lithium silicate. This set contains pre- and high-gloss polishers, which are optimally matched to each other in shape and grain. The 2-stage system can be used simply and without polishing paste on occlusal surfaces, cusps and fissures and demonstrates excellence with absolute high gloss.

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LUS04 LUSTER® for Zirconium reinforced Lithium Silicate

Intraoral Adjusting and Polishing Set

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

¹ 204=RA

LUS01 LUSTER® for Hybrid Ceramics Extraoral HP Laboratory Set

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The MEISINGER LUSTER® for Hybrid Ceramics polishing system has been specially developed for use on hybrid ceramic restorations. The system consists of two sets, which are adapted for intraoral and extraoral use. The polishing instruments provide pre-polishing and high gloss polishing for a coordinated two-stage system. Occlusal surfaces, cusps and fissures can be optimally processed without polishing pastes and are then impressive with their outstanding high gloss finish.

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LUS02 LUSTER® for Hybrid Ceramics Intraoral Adjusting and Polishing Set

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¹ 204=RA

GC05 Twist Polisher Lab Kit for Hybrid Ceramics

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The Twist Polisher Kit for Hybrid Ceramics is simple and efficient and was developed in collaboration with Greg Campbell, DDS, specifically to adjust and polish Hybrid Ceramics. The Lab kit contains an abrasive and diamond for sprue removal and refining and defining all surfaces. The sintered diamond S805-018 provides long service life and is excellent for enhancing occlusal anatomy. Twist Polishers specifically designed for Poly Ceramics bring out the luster of your restorations. The 9787-145 high shine wheel was specially developed to enhance the facial surfaces with stunning luster.

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GC06 Intraoral Twist Polisher Kit for Hybrid Ceramics

Intraoral Polishing Kit
according to Greg Campbell, DDS

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The Twist Polisher Kit for Hybrid Ceramics is simple and efficient and was developed in collaboration with Greg Campbell, DDS, to adjust and polish Hybrid Ceramics. The diamonds will gently adjust the occlusal surfaces and margins for proper fit and function. The polishing points and Twist Polishers will blend and polish all surfaces of the restoration.

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The Twist Polishing Kit has been especially developed for simple high-gloss polishing of all composite restorations. The flexible polishing discs make it possible to adapt to any surface structure and also render flawless results in areas that are difficult to reach, such as occlusal surfaces. Due to designs that are gentle on materials, the original surface structure of the restoration is maintained. In addition, the grades of diamond grit separated into two parts provide for short processing times and extremely high levels of durability.

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Diamond polishing system:

Especially for simple polishing of all composite restorations, also in regions difficult to access.

- Flexible polishing lamellae
- Specifically for occlusal surfaces
- Material-saving
- Short processing times
- Extremely high service life

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The MEISINGER Twist Finishers ensure a perfect working result and an extremely smooth surface on composites and other filling materials. Thanks to the unique cutting geometry, the instruments run extremely quietly and cut much finer than standard carbide finishers. The instruments are marked with colour rings for clear differentiation. The Twist Finishing Kit combines all available instruments in one assortment, with the lasered labels of the instruments on the Bur Block a sorting of the Twist Finisher is easily possible.

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MEISINGER LUSTER® Composite Finishing Kit was created to adjust and polish today's advanced composite materials. The contained MEISINGER Twist Finishing Burs were specifically designed for both hybrid and new generation composite materials. The helical design of the Twist Finishers ensures that the instrument stays in contact with the tooth surface to eliminate chatter, which results in an extremely smooth finish. The two-step polishing system allows for short processing times and long service life for the Twist Polishers, which are flexible enough to adapt to any surface structure while gentle enough to preserve that structure.

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These instruments will enable the achievement of the highest surface quality when finishing and polishing composite filling. In addition to the lingual processing of anterior teeth, the occlusal processing of molars and the processing of interproximal surfaces, the possibility of the selective processing of gaps in posterior teeth is also created. The various universally used polishers are an outstanding addition to this set. They also allow the perfect polishing of fissures and the finishing of the labial surfaces of anterior teeth.

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Flexible polishing discs for processing composites, ceramics and metals. This set contains 25 polishing discs of four different grits in Ø 10 mm and a carrier.

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Flexible polishing discs for processing composites, ceramics and metals. This set contains 25 polishing discs of four different grits in Ø 14 mm and a carrier.

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Three-step ceramic polishing system with knife-edge wheels and polishing points in coarse, fine and ultra-fine. These polishers are ideal for polishing ceramic in dental labs.

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LUS30 LUSTER® for Porcelain Intraoral Polishing Kit

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The LUSTER® Intraoral Porcelain Polishing Kit is specially developed for the polishing of porcelain ceramics. The instruments of this system are available in different grits, from coarse to extra fine, and allow rapid and systematic processing of all porcelain ceramics. This three-step polishing system produces high-shine results after a very short processing time. The use of polishing paste is not necessary.

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The last step in any indirect cementation appointment is the final polish. Dr. James Klim, DDS, developed the LUSTER® Intra Oral Polishing Kit to make this step as seamless and efficient as possible. These polishers come in three shapes and in four different grits.

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One kit for machining and polishing high-performance dental ceramics, such as aluminium oxide and zirconium oxide. This kit were specially designed for both the dental practitioner and the laboratory technician.

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
¹ 104=HP

LUS60 LUSTER® CAD/CAM

Porcelain Preparation and Polishing Kit

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
The kit includes all the preparation instruments for state-of-the-art preparation of CAD/CAM manufactured crowns and bridges. A highly polished result can even be achieved following corrections to the ceramic by using a 3-stage polisher.

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The Ceramics Polishing Kit is a three-tier system for polishing ceramics in dental laboratories. Diamond polishers in coarse, medium and fine grit are included, thus enabling a high-gloss finish in three steps.

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With the CAD/CAM Advanced-Kit you can carry out all stages necessary for preparation, adjustment and finishing of the CAD/CAM restoration quickly, efficiently and precisely.

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2662 Intra Oral Porcelain Adjustment Kit

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This set offers a wide selection of special instruments for intraoral processing, finishing and polishing of ceramic restorations. The 820F diamond grinder allows the safe and fast removal of excess material in the interdental region. The 890F fine-grained diamond was designed to improve the anatomy and interspaces in the molar region. Figure 833F offers the opportunity to efficiently finish both anterior teeth (palatal / lingual) and occlusal surfaces as well as molars. The set includes diamond ceramic polishers in two shapes and three different grain types for perfect polishing of the restorations.

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¹ 204=RA, 314=FG

KUNSTSTOFFE | ACRYLICS | ACRÍLICO

PMMA1 PMMA/PEEK Adjusting and Polishing Kit

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MEISINGER's PMMA / PEEK Adjusting and Polishing Kit includes all the necessary instruments for adjusting and polishing temporaries and partials. MEISINGER's fine cross-cut lab carbide HM75FX-040-HP for adjusting is included in this kit as well as MEISINGER's new polishers specifically designed for the polishing of PMMA and PEEK material. The cotton polishing brush completes the kit and provides a superior finish to the temporary/partial. The instruments are well organized in a high quality bur block.

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MEISINGER's Black Cobra Crown Temporization Kit includes all the instruments for adjusting and polishing temporaries. MEISINGER's Black Cobra fine lab carbides for adjusting acrylics are included in this kit. The unique Black Cobra carbides are notable for their increased longevity, a result of the low-friction coating that protects the carbide surface of the bur. Two diamond discs are included. The smaller disc is ideal for defining the embrasures, edges and corners. The three polishers complete the kit and provide a superior finish to the temporary.

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The Temporary Solution Kit is used for the simple processing and creation of provisionals in the dental practice.

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
The Acrylic Adjustment Kit correction set offers a comprehensive instrument kit for the quick and easy implementation of various corrections to synthetic materials in the practice. The innovative gear geometry of the hard metal burs included in the range of products thus already ensures the achievement of unusually smooth surfaces with the use of these instruments. The various synthetic material polishers are the result of years of development and are characterized by exceptionally high performance and durability.

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¹ 104=HP

VAL01 Flexible Acrylic Adjusting and Polishing Kit

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MEISINGER's Flexible Acrylic Adjustment and Polishing kit is the only kit you will ever need for adjusting and polishing Valplast®, DuraFlexTM and tcs® partials. The thin, flexible, and lightweight nature of these products requires rapid and efficient reduction without excess heat build-up. Lab carbides for both coarse and fine adjustments are included in addition to a fine polisher and cotton buffing wheel. The instruments are well organized.

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TITAN | TITANIUM | TITANIO

1321 Titanium Power Kit

developed with University RWTH Aachen

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UNCOMPROMISING for the processing of titanium.

The ST hard metal burrs are characterized by the super-coarse cutting geometry specially developed for titanium and an extremely rapid reduction in volume.

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Improve your endodontic access by improving visibility with Meisinger's Endo Access Bur Kit. This kit is a complete set of Meisinger Endo Access Burs designed for safe access into the root canal. Unlike traditional burs that obstruct your vision into the working area, these burs have long, thin necks allowing for an excellent, unobstructed view. These unique burs allow the practitioner to safely and easily remove dentinal triangles and flare the canal orifice for optimal straight line access into the root canal. Due to their very stable construction, the Endo Access Burs have optimal low-vibration concentricity and are designed to create less heat and debris than traditional ultrasonic tips. The color bands on the shank allow for quick identification and selection of the instruments.

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MEITRAC I is an endo-safety system for the precise gripping and fast, controlled removal of broken fragments such as root canal instruments, root filling posts, root canal posts and silver pins from the root canal. MEITRAC I is designed so that fragments with diameters of 0.15 to 0.50 mm can be safely removed.

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

The Periimplantitis Kit contains eight perfectly matched tungsten carbide finishing burs in egg and flame shapes for the ideal peri-implantitis treatment. Both shapes come in two sizes and each in standard and ultra-fine toothing. The shapes, sizes and toothings of the working parts are perfectly suited for the intraoral treatment of titanium and for different shapes of implant neck and shoulder. All finishers have an extra-long FG shank (total length 32 mm) which allows particularly the treatment of difficult-to-access deep areas. For the peri-implantitis treatment, the finishers are used clockwise and are guided around the exposed implant counterclockwise so the implant surface is smoothed.

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¹ 317=FGXXL

SONSTIGE | MISCELLANEOUS | OTRO

BC04 Black Cobra Complete Crown Removal Kit

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The BC04 Black Cobra Complete Crown Removal Kit is the only kit you will need for the removal of any dental restoration. Included in the kit are six Black Cobra diamonds for efficiently removing the hardest of crown materials such as Zirconia, as well as four special carbides for cutting through metal substructures or PFM crowns.

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Range of products with selected hard metal burrs of the MEISINGER Black Cobra Line for use in the dental laboratory. The special CABOCER® coating optimizes the performance of instruments and boasts extreme hardness, low wear and low friction values. The characteristic black color of the instruments allows you additional glare-free work.

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For removing residual plaque from subgingival areas. For treating teeth with two or three roots and severe bone recession which cannot be reached with a curette or scaler. Inaccessible areas should be extended with red ring diamonds. Deposits and calculus should be removed from the root surface with yellow ring diamonds. Best smoothing of root surfaces can be obtained with the white ring diamonds. Root surface treatment, when treating periodontitis, can be achieved minimally invasively utilizing this kit.

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MEISINGER's Fuzzies polishing wheels, when used with low speeds and light pressure, create a fine finish without overheating due to friction. The polishing wheels provide the user the ability to polish precious metals, composites and acrylics without damaging the surfaces.

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The Acryl/ Resin Universal Kit developed by Dr. James Klim, DDS, is especially suitable for the treatment and polishing of denture plastics, splints and aligners. With the included tungsten carbide instruments, it is easy to work out rails and dentures. The fiber fleece brushes provide a high-quality finish without overheating the material and damaging the surface. For a mark less surface, is the polisher 9573H included in the kit. All instruments are suitable for extraoral use and can be conveniently stored in the bur block.

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Create your own individual sets from the instruments you need for your daily work. The large selection of different MEISINGER drill stands offers you a high degree of flexibility when sorting your instruments. Not only are your dental instruments well organized in MEISINGER's high-quality drill stands, they can also be conveniently sterilized after use. The different drill cleaning brushes also ensure perfect cleaning.

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STAINLESS STEEL BUR BLOCK, NON-CORROSIVE, STERILIZABLE

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227

MEISINGER`S ELITE SERIES MIKROMOTOR | MEISINGER`S ELITE SERIES MICROMOTOR | MEISINGER`S ELITE SERIES MICROMOTOR

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MEISINGER`s micro-motor is equipped with a high-performance, extremely fatigue-free lithium ion battery to enable excellent performance.

Up to 5 hours of working time is possible on battery operation. Complete recharging only takes 3 hours. Despite its compact dimensions, the ultra-small micro-motor offers a performance superior to that of many stationary micro-motors: a maximum speed of 30,000 rpm and a robust 2.4 Ncm torque.

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These products are manufactured from colored anodized aluminum. When using chemical cleaning solutions, a slight discoloring may occur. Therefore, when using cleaning solution (recommended pH-value 5.5-8.5), it is advisable not to use those which have high contents of alkaline or acid.

The following sterilization method can be used without hesitation: Sterilizing with humid heat in the autoclave (e.g. at 134 °C / 273 °F). Please follow the operating instructions of the appliance according to the manufacturer and observe the recommendations regarding the duration of sterilization.

These Bur Blocks are only suitable for instruments with the specified shanks. Instruments with shorter shanks are not adequately secured and may fall out.

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The following application tips apply for aluminum products:

- Observe the regulations precisely of the disinfectant manufacturer: some are not suitable for aluminum!

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

Heißluft-Sterilisatoren: max. 180 °C / 356 °F

Rust preventing spray, sterilizable, for steel instruments

Purpose of application: The rust preventing spray prevents rusting during and after sterilization.

Applicable for all methods of sterilization.

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Bur cleaning brush, scrubbing
brush, brass

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Bur cleaning brush, scrubbing
brush, nylon

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Note: The cleaning brushes are not detachable. Therefore, brushes with residue that cannot be removed during cleaning should be discarded.

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- Content: 50 ml
- Compatible with a variety of metals and alloys
- Good lubricating effect
- Easy to remove with aqueous and solvent-based cleaners

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Processing (cleaning, disinfecting and sterilising)

of medical devices produced by Hager & Meisinger GmbH
Issued: June 2018

The medical devices produced and sold by Hager & Meisinger GmbH are reusable unless their label contains explicit information to the contrary. However, as a rule, it is the sole responsibility of the doctor/ expert using the devices to decide whether, depending on the respective case and the potential wear and tear of the products, he can re-use the products and how frequently he uses them. In case of doubt, it is always advisable to discard the products early and to replace them. The manufacturer Hager & Meisinger GmbH cannot guarantee the faultless function and performance of the products combined with a maximum degree of safety if the products are overused.

These reprocessing instructions apply in principle to all medical devices making up the product range supplied by Hager & Meisinger GmbH. Any particular features and/or exclusions that only concern individual items or groups of items are referred to separately. As to the general application and safety instructions concerning the use of the products, we would advise you to consult the application and safety instructions for the medical products of Meisinger, which are available separately (see also at www.meisinger.de).

General Principles

All products are to be cleaned, disinfected and sterilised prior to each application; this is required in particular for the first-time use after delivery of the unsterile instruments (cleaning and disinfecting after the removal of transport packaging; sterilisation after removing wrapping). An effective cleaning and disinfection is an indispensable requirement for an effective sterilisation of the instruments.

As you are responsible for the sterility of the products during use, please ensure

- that only sufficiently device and product specifically validated procedures will be used for cleaning, disinfection, and sterilisation
- that the used devices (disinfectant, sterilizer) are maintained and checked at regular intervals and
- that the validated parameters are adhered to during each cycle.

Please ensure to avoid a higher contamination of the complete sterilization tray during application by separate collection of contaminated instruments (without laying back into the sterilization tray). Pre-clean the contaminated instruments, then sort them back into the sterilization tray and clean, disinfect and sterilize the completely equipped sterilization tray.

Additionally, please pay attention to the legal provisions valid for your country as well as to the hygienic instructions of the doctor's practice or of the hospital. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA).

As the products are destined to be used for surgical, parodontological or endodontic procedures such as root canal debridement, they may penetrate the skin or the mucosa and come into contact with blood, internal tissues or organs (including wounds). Therefore, we recommend that they be assigned to risk group Critical B if used for their intended purpose. Attention: In the case of some instruments are additional or deviating procedures required (see chapter "Specific aspects").

Cleaning and disinfecting

Basic rules:

If possible, an automated procedure (WD (Washer-Disinfector)) should be used for cleaning and disinfection of the instruments. A manual procedure – even in case of application of an ultrasonic bath – should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of a manual procedure has to be considered. The pre-treatment step is to be performed in both cases.

Pre-treatment:

Abrasive impurities need to be removed from the products directly after use (within two hours maximum).

Procedure:

1. Disassemble the instruments as possible. Remove contaminated instruments of the sterilization tray (see chapter "Specific aspects").
2. Rinse the instruments at least 1 min under running water (temperature < 35 °C/95 °F) if applicable (see chapter "Specific aspects"): Rinse all lumens of the instruments at least three times at the beginning and at the end of the soaking time with a syringe (minimum volume 5-10 ml). Sway movable parts at least three times during pre-rinsing.
3. Soak the disassembled instruments for the given soaking time in the pre-cleaning solution¹ so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush¹ (at beginning of soaking, aids see chapter "Specific aspects") and subsequent ultrasonic treatment (after brushing, for the minimum soaking time, but not less than 5 min). Sway movable parts at least three times during pre-cleaning. If applicable (see chapter "Specific aspects"): Rinse all lumens of the instruments at least three times at the beginning and at the end of the soaking time with a syringe (minimum volume 5-10 ml).
4. If applicable (see chapter "Specific aspects"): Rinse all lumens of the instruments at least three times at the beginning and at the end of the soaking time with a syringe (minimum volume 5-10 ml).
5. Check the instruments on visible remnants. In case of still remaining remnants (e.g. bone or dentin particles) repeat steps 2 to 5, otherwise discard the instruments.

Pay attention to the instructions of the detergent manufacturer regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

¹ In case of application of a cleaning and disinfection detergent for this (e.g. in consequence of personnel's safety) please consider, that this should be aldehyde-free (otherwise fixation of blood impurities), possess a fundamentally approved efficiency (for example VAH/DGHH or FDA/EPA clearance or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material, and be compatible with the instruments (see chapter "material resistance"). Please consider, that a disinfectant used in the pre-treatment step serves only the personnel's

safety, but cannot replace the disinfection step later to be performed after cleaning.

2. Never clean products, bur blocks and sterilisation trays using metal brushes or steel wool.

Automated cleaning/disinfection (disinfector/ WD (Washer-Disinfector)):

Please consider the following points during selection of the WD:

- fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHH or FDA clearance)
- possibility for an approved program for thermal disinfection (A0 value ≥ 3000 or – in case of older devices – at least 5 min at 90 °C/194 °F; in case of chemical disinfection danger of remnants of the disinfectant on the instruments)
- fundamental suitability of the program for instruments as well as sufficient rinsing steps in the program
- post-rinsing only with sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water
- only use of filtered air (oil-free, low contamination with microorganisms and particles) for drying

If a WD is built in accordance with DIN EN ISO 15883 and regularly tested and maintained during its service life, it meets the above mentioned requirements with regard to water and air quality.

When choosing an appropriate cleaning and disinfecting agent you need to ensure

- fundamental suitability for the cleaning of instruments made of metallic or plastic material
- additional application – in case of non-application of a thermal disinfection – of a suitable disinfectant with approved efficiency (for example VAH/DGHH or FDA/EPA clearance or CE marking) compatible to the used cleaning detergent
- compatibility of the used detergents with the instruments (see chapter "material resistance,")

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing.

Procedure:

1. Disassemble the instruments as possible. Only pre-cleaned instruments can be sorted into the sterilization tray, e.g. bur blocks made of stainless steel (see chapter "Specific aspects").
2. Transfer the disassembled instruments in the WD (pay attention that the instruments have no contact). If applicable (see chapter "Specific aspects"): Connect the instruments to the rinsing port of the WD.
3. Start the program.
4. Disconnect (if applicable) and remove the instruments of the WD after end of the program.
5. Check and pack the instruments immediately after the removal (see chapters "check", "maintenance", and "packaging", if necessary after additional post-drying at a clean place).

The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the WD G 7836 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the pre-cleaning and cleaning detergent Neodisher mediclean forte (Dr. Weigert GmbH & Co. KG, Hamburg) considering to the specified procedure.

Manual cleaning and disinfection:

When choosing an appropriate cleaning and disinfecting agent you need to ensure

- fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- application of a disinfectant with approved efficiency (for example VAH/DGHH or FDA/EPA clearance or CE marking) compatible with the used cleaning detergent
- compatibility of the used detergents with the instruments (see chapter "material resistance,")

Combined cleaning/disinfection detergents should not be used. Only in case of extremely low contamination (no visible impurities) combined cleaning/disinfection could be used. Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

Procedure: Cleaning

1. Disassemble the instruments as possible. Only pre-cleaned instruments can be sorted into the sterilization tray (see specific dismantling instructions).
2. Soak the disassembled instruments for the given soaking time in the cleaning solution so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush¹ and subsequent with ultrasonic treatment (after brushing, for the minimum soaking time, but not less than 5 min). Sway movable parts at least three times during cleaning. If applicable (see specific dismantling instructions): Rinse all lumens of the instruments at least five times at the beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 5-10 ml).
3. Then, remove the instruments of the cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water. Sway movable parts at least three times during post-rinsing. If applicable (see specific dismantling instructions): Rinse all lumens of the instruments at least five times at the beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 5-10 ml).
4. Check the instruments (see chapters "check" and "maintenance").

¹ Never clean products, bur blocks and sterilisation trays using metal brushes or steel wool.

Disinfection

5. Soak the disassembled instruments for the given soaking time in the disinfectant solution so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Sway movable parts several times during disinfection. If applicable (see specific dismantling instructions): Rinse all lumens of the instruments at least five times at the

beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 5-10 ml).

6. Then, remove the instruments of the disinfectant solution and post-rinse them at least five times intensively (at least 1 min) with water. Sway movable parts at least three times during post-rinsing. If applicable (see specific dismantling instructions): Rinse all lumens of the instruments at least five times at the beginning and at the end of the soaking time by application of a single-use syringe (minimum volume 5-10 ml).
7. Dry and pack the instruments immediately after the removal (see chapter "packaging", if necessary after additional post-drying at a clean place).

The fundamental suitability of the instruments for an effective cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the cleaning detergent Gigazyme and the disinfectant Gigasept Instru AF (Schülke & Mayr GmbH, Nordstedt) considering the specified procedure.

Checking

After all products have been cleaned and/or cleaned/disinfected, check them for corrosion, damaged surfaces/ bare patches, broken/ chipped-off edges, deformations (e.g. bent rather than round) and impurities and eliminate damaged products (limited numbers for re-use see chapter on 'Re-use'). Products that are still contaminated need to be cleaned and disinfected once more.

Maintenance

- Re-assemble disassembled products (see specific instructions).
- Instrument oils must not be used.

Packaging

Please insert the cleaned and disinfected products in the dedicated bur block/sterilisation tray. Please pack the instruments or the sterilization trays single-use sterilization packaging (single or double packaging) and/or sterilization containers, which fulfill the following requirements (material/process):

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to at least 138 °C (280 °F), sufficient steam permeability)
- sufficient protection of the instruments as well as of the sterilization packaging to mechanical damage
- regular maintenance according to the instructions of the manufacturer (sterilization container)
- Individual packaging: the packaging must be sufficiently large to ensure that the sealing is tension-free.

Sterilisation

We only recommend the use of the sterilisation procedures listed below.

Steam sterilization:

- fractionated vacuum/dynamic air removal procedure¹ or gravity displacement procedure² (with sufficient product drying)³
- steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- maximum sterilization temperature 134 °C (273 °F; plus tolerance according to EN ISO 17665)
- sterilization time (exposure time at the sterilization temperature):

Area	fractionated vacuum/dynamic air removal	gravity displacement
USA	at least 4 min at 132 °C (270 °F), drying time at least 20-30 min ⁴	at least 4 min at 132 °C (270 °F), drying time at least 15-30 min ⁴ at least 30 min at 121 °C (250 °F), drying time at least 15-30 min ⁴
other countries	at least 3 min ⁴ at 132 °C (270 °F)/134 °C (273 °F) at least 20 min at 121 °C (250 °F)	at least 5 min at 132 °C (270 °F)/134 °C (273 °F) at least 20 min at 121 °C (250 °F)

¹ at least three vacuum steps

² The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure.

³ The effectively required drying time depends directly on parameters in sole responsibility of the user (load configuration and density, sterilizer conditions, ...) and by this is to be determined by the user. Nevertheless, drying times less than 20 min must not be applied.

⁴ respectively 18 min (inactivation of prions, not relevant for USA)

The fundamental suitability of the instruments for an effective steam sterilization was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the steam sterilizer SysTec V-150 (SysTec GmbH Labor-Systemtechnik, Wetzlar) and the fractionated vacuum/dynamic air removal procedure as well as the gravity displacement procedure. For this, typical conditions in clinic and doctor's practice as well as the specified procedure were considered. The flash/immediate use sterilization procedure must not be used. Do not use dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization.

Storage

Prior to the first use of the device, the product should be stored in its original packaging at room temperature in dust- and humidity-free conditions. Subsequently, the products should be stored in appropriate hygienically maintained containers (protected from dust, humidity and recontamination). After sterilisation, the products need to be stored in sterilisation wrapping in a dry and dust-free place. Please note the shelf-life resulting from the validation of the sterilisation wrapping.

Material resistance

When choosing the cleaning and disinfecting agents ensure that they do not contain the following ingredients:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- strong lyes (maximum admitted pH-value 11, neutral/enzymatic, weak alkaline or alkaline cleaner recommended)
- organic solvents (for example: acetone, ether, alcohol, benzene)
- oxidizing agents (for example: peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons

Please consider during selection of the detergents in addition, that corrosion inhibitors, neutralizing agents, and/or rinse aids may cause potential critical remnants on the instruments. Please do not clean any instruments and sterilization trays by use of metal brushes or steel wool. Please do not expose any instruments and sterilization trays to temperatures higher than 138 °C (280 °F)!

Re-use

The instruments can be reused – in case of adequate care and if they are undamaged and clean. For bur blocks – in case of adequate care and if they are undamaged and clean – a reuse up to 100 times is possible. The user is responsible for each further use as well as for the use of damaged and dirty instruments (no liability in case of disregard).

Specific aspects

Benex-Control (Extractor):

- Disassemble instrument completely (including the Teflon plate) prior to cleaning and disinfecting
- Reassemble instrument prior to wrapping and sterilisation; ensure that the ball does not catch in a cavity of the biteplate (semi-tilted position)

Benex-Control (Traction cable):

- Bend cable during pre-cleaning and manual cleaning and disinfection several times in all directions

Diamond products and ceramic grinding tools:

- Use particular care during the cleaning of the grinding surfaces and ensure that all residues are removed

Distractors:

- Move instrument up and down several times during pre-cleaning and manual cleaning and disinfecting
- Mechanical cleaning/disinfection as well as sterilisation only in specific positions (VCD05: half-open state, VCD15: nearly closed state)

Products equipped with a conduit to supply a cooling medium and other products with lumens (canals, drillings etc.):

- Actively rinse the lumen during pre-cleaning and manual cleaning and disinfection process
- Products containing lumens without feed-through channels must not be re-used

Bur blocks / instrument trays:

- Bur blocks made of stainless steel (Article no. BS440, BS470, BS540, BS570, BSS37 und BBS55) can be cleaned / disinfected loaded with instruments in an automated WD (washer / disinfectant).
- Bur blocks made of other materials, e.g. aluminum must not be cleaned and disinfected whilst they are loaded with instruments. Remove attachments prior to cleaning and disinfecting and clean and disinfect in disassembled state.

Universal Torque Ratchet:

- Disassemble instrument completely prior to cleaning and disinfecting
- Move joint several times to and fro during pre-cleaning and manual cleaning and disinfecting
- Reassemble instrument prior to wrapping and sterilisation.

Instruments made of tool steel (tungsten-vanadium):

- Instruments made of tool steel are not suitable for automated cleaning/disinfection.
- Please consider for manual cleaning and disinfection also the chapter material resistance. For manual cleaning a neutral/ enzymatic detergent is recommended.
- Without appropriate pre-treatment the instruments are not suitable for sterilization.

For additional information on the procedure, please refer to the notes on „Manual processing of Hager & Meisinger GmbH instruments made of tool steel“ in the download area under www.meisinger.de

The above instructions have been validated by the manufacturer of the medical devices who found them to be SUITABLE for preparing a medical device for re-use. It is up to the person in charge of the reprocessing to ensure that, based on the use of the correct equipment, material and personnel in the reprocessing facility, the actual reconditioning process produces the desired results. Normally, this requires the validation and routine monitoring of the procedure. Equally, each deviation from the instructions provided should be carefully checked for effectiveness and potential adverse consequences by the person in charge of reprocessing.

General application and safety instructions for the medical device produced by MEISINGER

Issued: October 2018

- MEISINGER medical products (dental, implantology, Bone Management, maxillary surgery, general surgery, chiroprody devices) must only be used by dentists, doctors and/or the respective experts who, due to their training and experience, are intensely familiar with the use of these products and who have the corresponding expertise in the respective specialist fields. The use of surgical products requires relevant expertise and experience in dental implantology, maxillary surgery and/or other surgical fields including diagnosis, preoperative planning and surgical techniques.
- It is the sole responsibility of the doctor in charge who, depending on the respective situation (indication), decides on the actual use of the MEISINGER products for each individual case.
- Courses providing further training on the possible uses of MEISINGER products are organised on a regular basis. You will find relevant information about these events as well as further information on MEISINGER products on the internet at www.meisinger.de.
- All MEISINGER products have been developed for specific applications. Therefore, inappropriate use can result in the premature wear and tear of the products and put patients and users at risk.

Application

- In order to avoid damaging the instruments, they must be removed from the blister pack individually.
- It is essential to only use turbines as well as hand and angle pieces that are technically and hygienically faultless, maintained and cleaned.
- The instruments must be rotating when applied on material. They should not be placed on material and then brought to rotation.
- Rotating instruments need to be clamped as far down as possible with their speed set before applying them on the object. are used with the rotary instruments.
- Using the instruments for canting or leveraging should be avoided as it increases the risk of breakage.
- Depending on the application, it is recommended to use protective goggles while using the instruments. Users of diamond disks should use a disk protection device.
- Inappropriate use of the products leads to badly executed work and increased risk.
- When working with dry materials, it is recommended to use a suction cleaning device.
- In particular, users of hand tools should take care to use them gently and with consideration.
- The user must at all times avoid touching the instruments and parts without protection (protective gloves should be worn).
- Thermal bone damage caused by rotating and oscillating tools (e.g. pilot burr, conical burr, expansion burr) must at all times be avoided (user training, working at low speed and with sufficient cooling).
- During intraoral application attention has to be made to the fact that the products are protected against aspiration or falling on the floor.

Use of pressure

- Users of the instruments should at all times avoid applying excessive pressure. This can damage the working part of the instruments and cause the cutting edges to break off. At the same time, it generates excessive heat.
- The use of excessive pressure when using grinding tools can cause the abrasive particles to break off or the instrument to become clogged and lead to heat generation.
- During polishing, excess pressure can lead to heat generation.
- Due to overheating, excess pressure can damage the dental pulp or, due to broken off cutting edges, it can result in undesired rough surfaces. In such cases, even instrument breakage cannot be excluded.

Cooling

- In order to avoid excessive heat generation during preparation, a sterilewater/sodium chloride solution supplied via a permanent external feeding device should be used to ensure sufficient cooling during use of the instruments.
- When using FG instruments that are more than 22 mm long or whose head diameter exceeds 2 mm, additional external cooling is required.
- Insufficient cooling will lead to irreversible damage to the bone and/or the adjacent tissue.

Storage, disinfection, cleaning and sterilization

- Unless there is explicit information to the contrary, all MEISINGER products are supplied in non-sterile packaging and, depending on the application, they need to be sterilised prior to use. Prior to their first use on the patient and immediately after each use, all products need to be disinfected and sterilised. Inappropriate cleaning and sterilising of the instruments can result in the patient being infected with harmful bacteria.
- You will find detailed instructions for the disinfecting, cleaning and sterilising of products in the Instructions for the processing of medical devices produced by Hager & Meisinger GmbH. We would be happy to provide you with these instructions at your request. They are also available on the internet at www.meisinger.de.
- The products should be stored in appropriate, hygienically maintained containers. The same applies to sterilised instruments. The stored products must be protected from dust, humidity and recontamination. Instructions as to maximum storage duration must be adhered to.

Speed recommendations for rotary instruments

- Following the instrument-specific speed recommendations produces the best results.
- Exceeding the maximum admissible speed (rpm) when using long and pointed instruments tends to produce vibrations that can lead to the destruction of the instrument.
- When using working parts with diameters exceeding the thickness of the shaft, excessive speed can release great centrifugal forces that may cause the shaft to bend and/or the instrument to break. Therefore, the maximum admissible rpm must never be exceeded.
- Please consult the manufacturer's information (see catalogue or www.meisinger.de) for the recommended and the maximum admissible speed ranges. Non-compliance with the maximum admissible speed puts safety at risk.
- Generally, the following rules apply:
 - The larger the working part of an instrument the lower the speed
 - The larger the working part of an instrument the greater the force of pressure
- Instruments with a maximum speed of 300,000 rpm are suitable for micro-motor hand pieces and turbines with stable ball bearings. Not to be recommended for turbines with air bearings.

- Instruments with a maximum speed of 30,000 to 160,000 rpm are suitable for micro-motor hand pieces or technical hand pieces up to the indicated rpm. Not to be recommended for turbines.
- Surgical instruments: suitable for geared down micro-motor hand and angle pieces 10:1 with stable ball bearings. Speed 600 to 800 rpm with physical and, possibly, sterile external cooling or internal cooling when using the respective hand piece.

Discarding worn instruments and parts

- MEISINGER products can principally be reused several times – unless specifically indicated and labelled otherwise. Rotating instruments are subject to wear. The option of and accountability for multiple use of a product and the frequency of application is solely the decision and own responsibility of the treating clinician based on the application in each case and the possible wear of the products. If in doubt, the products should always be sorted out early and replaced.
- Broken off cutting edges of instruments cause vibrations and great forces of pressure, which, in turn, leads to broken preparation corners and rough surfaces.
- Bare patches on diamond instruments indicate a lack of abrasive particles and can be a sign of blunt instruments. This leads to excessive temperatures during instrument use.
- Instruments that are bent and/or do not run true should be discarded forthwith.
- With the reuse of disposable products the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed.

Additional instructions for the use of implants

- As a rule, implants are meant to be used only once.
- Implants must only be inserted by dentists and doctors who are intensely familiar with dental implantology including diagnosis, pre-operative planning, surgical techniques and prosthetic care.
- The unprotected touching of the implant by the user must be avoided (protective gloves must be worn).
- Systemic errors that are due to a mix-up of tools and implants are to be avoided. Therefore, paying attention to colour coding and labelling is essential.
- A lack of primary stability can be due to the faulty preparation of the implant bed as well as to deficient bone quality. In such cases, the mobile implant should be removed and replaced with a larger one.
- The best way to avoid complications is to only use implants when all the conditions for a complication-free treatment guaranteeing long-term success have been met. If, from the start, there are reservations regarding the chances of success of a desired implant-based dental prosthesis, it may be best to choose another prosthetic solution that does not require an implant.

Additional instructions regarding the use of trepans

- When using trepans, you have to proceed with particular care. For example, it is advisable not to exceed the recommended rpm speed ranges.
- In order to prepare for the actual use of a trepan, it should be set to produce counter-clockwise rotations creating a groove in the bone. Afterwards the trepan can be inserted into this groove and, using clockwise rotations, it can be moved further down.
- Carrying out a prior X-Ray is essential to establish the maximum possible drilling depth and to maintain the necessary distance, for example, to the mandibular nerve. As an additional safety measure to spare the nerve, the axial direction of the trepan countersink attachment, based on the sagittal level of the ascending branch, must be milled laterally at an angle of approx. 15-20°.

Additional instructions for the use of oscillating instruments

- Prior to each use, the instruments should be tested for operational safety by applying slight pressure with the index finger or the thumb (risk of instrument fracture through continuous use or due to the instrument having been dropped). Prior to use, the instrument must be firmly screwed down because, otherwise, reduced performance may be the outcome.
- Complying with instrument-specific performance settings (frequency, amplitude of vibration) produces the best results. Non-compliance with recommended settings may lead to premature breakage of the instruments.
- The instruments are generally suitable for ultrasonic transmitters that provide frequencies between 24 kHz and 36 kHz for surgical applications or in the field of 5 to 7 kHz for prophylactic applications with amplitudes of vibration between 60 and 300 µm.
- Only such drive mechanisms should be used that are designed in a way that sufficient cooling can be provided to the prepared site.

Additional information for the implementation of the SINGLES (supplied sterile single use product)

- The SINGLES are delivered sterile and are for single use only
- The instruments are sterile as long as the bag remains unopened and undamaged. If the inner packing has been damaged in any way, bacteria contamination may have taken place which means that the instruments cannot be used nor can an unauthorized sterilization be carried out. Do not use the products after the specified date of expiry has elapsed (symbol hourglass).
- To remove the instruments, tear open along perforation on side of bag, on no account should the instruments be pressed out of the sterile packing.

Further comments:

- Due to statutory regulations, returned goods can, on principle, only be accepted if the complete batch number is provided. This number can be found on the product packaging.

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Rotating instruments by MEISINGER cover the entire spectrum of applications in dentistry. The respective speed is a theoretical value given on the basis of the geometry of the instruments and the principal application. You can find the specific maximum and recommended speeds in the product-related information.

- Exceeding the maximum permissible speed can lead to an impairment in the performance of the products and/or to a risk for the patient or the user.
- In particular, the long and sharp instruments tend to vibrate if the maximum permitted rotational speed is exceeded and this can lead to destroying the instruments.
- In the case of workpiece diameters which exceed the shaft strength, excessive speeds can lead to pronounced centrifugal forces which can lead to distortion of the shaft and/ or breaking of the instrument. For this reason, the admissible maximum number of revolutions must never be exceeded.
- The recommended speeds and maximum permissible speeds are given in the manufacturer's instructions (see brochure and product-related packaging).
- Please observe adapting the rotation speed in relation to the prevailing indication and the material as well as the diameter of the working part and the type of tooth connection. The selection of a working speed for a concrete application depends specifically on the instrument drive used as well as the maximum application forces and lies within the discretion of the user.

As a general rule:

- the larger the working part of an instrument, the lower the speed should be set.
- the larger the working part of an instrument, the greater the application force.
- the rougher the tooth connection of the article, the lower the selected speed should be.
- Instruments with a maximum speed 300,000 rpm: suitable for micro-motor handpieces and turbines with stable ball-bearings. Not recommended for air-cushioned turbines.
- Instruments with a maximum speed 30,000–160,000 rpm: suitable for micro-motor handpieces and technical handpieces up to the given speed. Not recommended for turbines.
- Surgical instruments: suitable for geared-down micro-motor handpieces and angled handpieces 10:1 with stable ball bearings, if applicable, with external cooling or internal cooling when using the appropriate handpiece.

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300 - 500 min⁻¹

Special notes**Diamond Instruments Practice:**

On hard materials such as dental enamel, it is recommended to use a finer grain at a higher speed for precision machining of the surface. The final decision as to graining and speed lies with the user.

Diamond Instruments Laboratory:

For delicate work, the speed should be selected according to the surface quality and the intended result of the work, taking into account the experience of the user. Please note that diamond instruments develop their best cutting quality at higher speeds due to the self-sharpening effect. The final decision as to graining and speed lies with the user.

Ceramic grinders

For hard materials (e.g. fired ceramics, precious metal (Co-Cr-Mo alloys) the use of soft grinding tools at high speeds is recommended to achieve abrasion but not „smudging“. For soft materials, machining with hard-bonded grinding tools at low speeds is recommended.

Surgical instruments

The following reference values for rotation speeds apply to surgery in general:

Handpiece

- recommended: 6,000 - 10,000 rpm
- maximum: 40,000 - 50,000 rpm

Angled handpiece:

- recommended: 6,000 - 10,000 rpm
- maximum: 40,000 - 50,000 rpm

FG:

- recommended: 80,000 rpm
- maximum: 100,000 - 120,000 rpm

231DC:

- recommended: 800 rpm
- maximum: 1,000 rpm

Trepans

Special care should be exercised when using trepans. Specifically, the recommended speeds are not to be exceeded. When preparing the use of a trepan, this should be inserted counter-clockwise to create a groove in the bone. Then the trepan can be inserted deeper into groove with clockwise rotation.

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All rights reserved.
Reproduction, also by extract and reproducing photos, are only permitted with written authorization of Hager & Meisinger GmbH, Neuss.
All offers, orders, and deliveries are subject to Meisinger's "General Sales and Delivery Terms."
We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

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