

NEW improved delivery system
NEW lower cost mixing tip refills.¹



Same great material.

Less waste. Improved mixing. Easy to use.

KEY SYSTEM CHANGES:

New Innovative Mixing Tips

- RED mixing tips for Low Viscosity material
- BLUE mixing tips for High Viscosity material
- New lower cost mixing tip refills¹
- New mixing tip item numbers
- **ONLY compatible with new cartridge system**

New Cartridges










- New color coded cartridge
- RED - Low Viscosity
- BLUE - High Viscosity
- Use with your current dispenser

Material

- No changes
- Same great material
- Same re-order numbers




Contact your Dentsply Sirona Representative for more information

Order Numbers for Old Mixing Tips		Order Numbers for New Mixing Tips	
Item #	Description	Item #	Description
<i>for Low Viscosity Material (Wash Material) XLV, LV</i>			
678212	Teal Mixing Tips 	778215M	RED Mixing Tips 
619451	Yellow Mixing Tips 		
626382	Clear Intraoral Mix Tips 	626383	RED Intraoral Tips 
<i>for High Viscosity Material (Tray Material) Medium, Heavy, Rigid, & Putty</i>			
678212	Teal Mixing Tips 	778213M	BLUE Mixing Tips 
678950	Dynamic Mixing Tips (Blue) 	778951M	Dynamic Mixing Tips 

New Cartridge Delivery System

- Can eliminate **up to 64%** of the waste you're used to²
- Bleed once for the life of the cartridge³
- **Up to 10 more** single-unit crown applications per cartridge²



1. When comparing current Dentsply Sirona mixing tip refills to new mixing tip refills
2. When comparing Dentsply Sirona low viscosity impression material in the new red mixing tip versus the old teal mixing tip.
3. Only with the 50mL System

ML780109F (7-12-21)