


safco Surestitch™ Sutures

Dental Supply




Plain gut (catgut)

Made of first quality raw material, catgut has extremely high tensile strength. Every strand is precision polished to a uniform diameter, permitting smooth and secure knotting.




Black braided silk

Noncapillary silk has excellent handling and tying characteristics. Modern braiding technique provides a uniform smooth surface and a greater tensile strength.



Chromic gut (catgut)

Produced with the same precision and quality materials as plain gut, plus treated with chromic acid salts to decrease tissue reactivity and provide a slower absorption rate.



Nylon black monofilament

A non-absorbable polyamide suture with characteristically high tensile strength. The uniform smooth surface permits easy passage through tissue.



Needle Type	Material	Safco code	Size	Length	Pkg size
C-3 P-3 	Chromic gut	#102-5865	5-0	18"	12/box
C-6 FS-2 	Black braided silk	#102-5833	4-0	18"	12/box
	Black braided silk	#102-5847	3-0	18"	12/box
	Chromic gut	#102-5801	4-0	18"	12/box
	Chromic gut	#102-5807	4-0	27"	12/box
	Chromic gut	#102-5811	3-0	18"	12/box
	Chromic gut	#102-5815	3-0	27"	12/box
	Nylon black monofilament	#102-5885	4-0	18"	12/box
	Plain gut	#102-5822	4-0	18"	12/box
	Plain gut	#102-5879	4-0	27"	12/box
	Plain gut	#102-5870	3-0	18"	12/box
Plain gut	#102-5829	3-0	27"	12/box	
C-7 FS-1 	Plain gut	#102-5875	3-0	27"	12/box
C-31 X-1 	Black braided silk	#102-5859	3-0	18"	12/box
	Plain gut	#102-5861	3-0	18"	12/box

■ Plain gut
 ■ Chromic gut
 ■ Black braided silk
 ■ Nylon black monofilament