

Finishing Strips

Abrasive and serrated strips designed for interproximal finishing on enamel or restorative materials, as well as breaking contact or creating interproximal reduction.

Most Popular

SmartStrip™ Finishing Strips


IPR • Cement Removal • Finishing

KaVo Kerr SmartStrip Finishing Strips are designed to make cleanup quick, safe, and simple after cementation of crowns, veneers, bridges, inlays or onlays. SmartStrips are also great tools for precise interproximal reduction with sequential gauges and accompanying technique guides. The unique handle provides a comfortable grip for optimal tactile control, giving easy access to interproximal spaces.



- Easy interproximal reduction
- One-handed operation
- Easy access to tight anterior and posterior spaces
- Large grinding surface and improved working view*
- Autoclavable, multiple-use device

Serrated - 10 per pack



Grit THICKNESS* (mm)	Serrated 0.05
MFG. NO.	SS-SS

Single-Sided - 10 per pack



Grit THICKNESS* (mm)	Coarse 0.18	Medium 0.14	Fine 0.10	Superfine 0.08
MFG. NO.	SS-C	SS-M	SS-F	SS-SF
Assortment	SS-ASST (2 of each and (2) SS-SS)			

Double-Sided- 10 per pack



Grit THICKNESS* (mm)	Coarse 0.35	Medium 0.22	Fine 0.20	Superfine 0.18
MFG. NO.	SSDS-C	SSDS-M	SSDS-F	SSDS-SF
Assortment	SSDS-ASST (2 of each and (2) SS-SS)			

Curved Single-Sided- 10 per pack



Grit THICKNESS* (mm)	Coarse 0.18	Medium 0.14	Fine 0.10	Superfine 0.08
MFG. NO.	SSC-C	SSC-M	SSC-F	SSC-SF
Assortment	SSC-ASST (2 of each and (2) SS-SS)			


* thickness is approximate

■ 120-135 Micron ■ 100-120 Micron ■ 45 Micron ■ 35 Micron


Most Popular

NTI® Serrated Diamond Finishing Strips

FS3-X Wide Serrated Series - 3.75 mm wide by 147 mm long



The NTI® Serrated Diamond Finishing Strip from KaVo Kerr Dental features an ultra-thin strip and a 40mm diamond-free serrated edge for breaking interproximal contacts. The diamond coated sections are for the finishing of enamel and all restorative materials. Autoclavable.




- Medium grit (37-44 Micron), 0.13 mm thick, 10 pack FS3-MX
- Fine grit (30-40), 0.10 mm thick, 10 pack FS3-FX
- Superfine grit (10-20 Micron), 0.08 mm thick, 10 pack FS3-SFX

Most Popular

Perforated Diamond Finishing Strips

2.5 - 4.0 mm wide by 150 mm long

- Diamond-free zone for easier fit and preparation
- Stainless steel strip resistant to breaking and stretching
- Autoclavable - reusable for greater value



- 2.5 mm, Medium grit (37-44 Micron), 0.11 mm thick, 10 pack FSP2-M
- 2.5 mm, Fine grit (30-40), 0.10 mm thick, 10 pack FSP2-F
- 2.5 mm, Superfine grit (10-20 Micron), 0.08 mm thick, 10 pack FSP2-SF
- 4.0 mm, Medium grit (37-44 Micron), 0.11 mm thick, 10 pack FSP4-M
- 4.0 mm, Fine grit (30-40), 0.10 mm thick, 10 pack FSP4-F
- 4.0 mm, Superfine grit (10-20 Micron), 0.08 mm thick, 10 pack FSP4-SF

Assortment Pack FSP42-SFMX:
 NTI® Serrated: FS3-MX (2 ea) / FS3-FX (1 ea) / FS3-SFX (1 ea)
 Axis 2.5mm: FSP2-M (1 ea) / FSP2-F (1 ea) / FSP2-SF (1 ea)
 Axis 4.0mm: FSP4-M (1 ea) / FSP4-F (1 ea) / FSP4-SF (1 ea)

NTI® Narrow Diamond Finishing Strips


FS2 Series - 2.5 mm wide by 147 mm long



- Medium grit (37-44 Micron), 0.13 mm thick, 10 pack FS2-M
- Fine grit (30-40 Micron), 0.10 mm thick, 10 pack FS2-F
- Superfine grit (10-20 Micron), 0.08 mm thick, 10 pack FS2-SF

NTI® Wide Diamond Finishing Strips

FS3 Series - 3.75 mm wide by 147 long



- Medium grit (37-44 Micron), 0.13 mm thick, 10 pack FS3-M
- Fine grit (30-40 Micron), 0.10 mm thick, 10 pack FS3-F
- Superfine grit (10-20 Micron), 0.08 mm thick, 10 pack FS3-SF

Wide Finishing Strip and Assortment Pack

- Medium grit (37-44 Micron), 6 mm wide, 10 pack FS6-M
- Assortment Pack contains 2 FS3-MX, 1 each of FS3-FX, FS3-SFX, FS32-SFMX, FS3-M, FS3-F, FS3-SF, FS2-M, FS2-F, FS2-SF

Diamond Finishing Strip Tube

165mm x 19mm..... AT600-1B



DOUBLE-SIDED ORTHO-STRIPS

4 per pack color coded by grit

REF. NO.	OS-15	OS-25	OS-40	OS-60
THICK mm	0.10	0.12	0.13	0.30
Grit	15	25	40	60
Grit Micron	15	25	40	60
Assortment	OS-SET-4 : 1 of each OS-15,25,40,60 (4 Pack)			
Package	4 per pack			



SINGLE-SIDED ORTHO-STRIPS

1 per pack color coded by grit

REF. NO.	OS-40L	OS-60L	OS-40R	OS-60R
THICK mm	0.12	0.13	0.12	0.30
Grit	40	60	40	60
Grit Micron	40	60	40	60
Assortment	OS-SET-SS6 : 1 of each OS-15,25,40, L & R (6 Pack)			
Package	Individual			




Intensiv Ortho-Strip assortment pack; 1 of each size .10mm - .30mm. MFG. NO.: OS-SET-4

Intensiv Ortho-Strip assortment pack; 1 of each side, size .10mm - .12mm. (6 pack) MFG. NO.: OS-SET-SS6

Intensiv™ Proxoshape

Follow and finish your procedure with Intensiv Proxoshape Diamond Files. Use these flexible Intensiv Proxoshape files to contour and polish tooth enamel. These files are also perfect for removal of amalgam or composite overhangs. Intensiv Proxoshape Diamond Files are available in numerous lengths and grits allowing efficient interproximal contouring that cannot be achieved with diamond burs, discs or strips.



The use of Intensiv Ortho-Strips and ProxoShapes require a specialhandpiece attachment sold by KaVo. Please contact your KaVo or distributor rep for more information.




3 per pack color coded by grit

REF. NO.	PS-0	PS-1	PS-2	PS-3	PS-1L	PS-2L	PS-3L	PS-2S	PS-3S	PS-3G
THICK mm	0.45	0.35	0.12	0.10	0.35	0.12	0.10	0.12	0.10	0.10
Grit	125	90	40	15	90	40	15	40	15	15
Grit Micron	125	90	40	15	90	40	15	40	15	15
Assortment	PS-100-4 : 1 of each plus File Injector & File Cleaner					PS-115-10 : 1 of Each				
Package						2 Per Pack				













Intensiv Proxoshape assortment pack; 1 of each size. (10 pack) MFG. NO.: PS-115-10

KaVo Kerr Rotary Polishers are the choice of dental professionals worldwide and are offered in a variety of shapes, sizes and shanks for adjusting and polishing on ceramics (CeraGlaze®), composites (PDQ™), acrylics, alloys, and gold. It is recommended that Polishers be used in sequence with recommended speeds (see below) and lighter pressure to achieve the highest polish possible. Polishing pastes are also available for the final polish on ceramics or composite.

Shank Type

Code	Shank Type	Image	Diameter (mm)	Length (mm)
FG	Friction Grip		1.6	19
RA	Right Angle Latch		2.35	22
HP	Handpiece		2.35	44.5

Polisher Application Guide

	Most Popular										
	OptiDisc	ProGloss™	Cera Glaze®	UltraGloss™	PDQ2™ Polishers	PDQ™ brush	DiaGloss®	Brown/ Green Silicone	Goldstar	White/ Pink/ Gray	Zir-Cut™
											
Application	• adjusting • pre-polish • high shine	• high shine	• pre-polish • polish • high shine	• adjusting • pre-polish • high shine	• high shine	• high shine	• pre-polish • high shine	• pre-polish • polish	• adjusting • pre-polish • high shine	• adjusting • pre-polish • polish	• pre-polish • polish • high shine
Shapes	disc	flame, cup, square edge disc	disc, knife edge/square edge disc, flame, cup	disc	flame, cup, square edge disc, small point	cup, point	disc, flame, cup, square edge disc	disc, flame, cup, knife edge/square edge disc	disc, flame, knife edge disc	disc, flame, cup, knife edge disc	disc, point, cup
Available shank	RA	RA	RA, FG, HP	RA	RA	RA	RA, FG	RA, FG, HP	RA, HP	RA, HP	RA
Recommended RPM	10,000 to 12,000	5,000	green - 16K* blue - 12K* yellow - 6K*	10,000 to 12,000	5,000	2,000	10,000 to 15,000	10,000 to 15,000	green - 20K yellow - 15K pink - 10K	white - 15K pink - 10K green - 5K	7,000 to 10,000
Amalgam		★★			★★			★★★			
Composite	★★★	★★★★		★★★	★★★	★★★	★★★	★			
Monolithic Zirconia			★★								★★★
Alloys non-precious and semi-precious			★★★					★★★	★★		
Gold gold and other precious alloys			★★★					★	★★★		
Porcelain			★★★							★★	

★ Good ★★ Better ★★★ Best
















*All CeraGlaze® FG shanks 10,000-12,000 RPM and RA shank 5,000-15,000 RPM. CeraGlaze®, and DiaGloss® are registered trademarks of NTI®.


Polishers

KaVo Kerr Rotary Polishers are the choice of dental professionals worldwide and are offered in a variety of shapes, sizes and shanks for adjusting and polishing on ceramics (CeraGlaze®), composites (ProGloss™), acrylics, alloys, and gold.






















DiaLite® to KaVo Kerr/NTI® CeraGlaze® Conversions

DiaLite®

 W16DG Intra-oral RA/Latch Bulk Reduction	 W16DM Intra-oral RA/Latch Prepolish	 W16D Intra-oral RA/Latch High Shine
 W18DG Intra-oral RA/Latch Bulk Reduction	 W18DM Intra-oral RA/Latch Prepolish	 W18D Intra-oral RA/Latch High Shine
 W17DG Intra-oral RA/Latch Bulk Reduction	 W17DM Intra-oral RA/Latch Prepolish	 W17D Intra-oral RA/Latch High Shine
 H2DG Extra-oral Handpiece / Bulk Reduction	 H2DM Extra-oral Handpiece / Prepolish	 H2D Extra-oral Handpiece High Shine
 R17DG Extra-oral Handpiece / Bulk Reduction	 R17DM Extra-oral Handpiece / Prepolish	 R17D Extra-oral Handpiece / High Shine

 L26DBC Knife Edge Handpiece Bulk Reduction
 L26DRM Knife Edge Handpiece Prepolish
 L26DGSF Knife Edge Handpiece High Shine

KaVo Kerr/NTI® CeraGlaze®

 P332-3 Intra-oral RA/Latch Prepolish	 P3032-3 Intra-oral RA/Latch Refining	 P30032-3 Intra-oral RA/Latch High Shine
 P336-3 Intra-oral RA/Latch Prepolish	 P3036-3 Intra-oral RA/Latch Refining	 P30036-3 Intra-oral RA/Latch High Shine
 P335-3 Intra-oral RA/Latch Prepolish	 P3035-3 Intra-oral RA/Latch Refining	 P30035-3 Intra-oral RA/Latch High Shine
 P344 Extra-oral Handpiece Prepolish	 P3044 Extra-oral Handpiece Refining	 P30044 Extra-oral Handpiece High Shine
 P343 Knife Edge & Wheel Extra-oral Handpiece / Prepolish	 P3042 Knife Edge & Wheel / Extra-oral Handpiece / Refining	 P3043 Knife Edge & Wheel / Extra-oral Handpiece / High Shine
 P310 Knife Edge Premounted Prepolish	 P301 Wheel Premounted Prepolish	
 P3010 Knife Edge Premounted Refining	 P3001 Wheel Premounted Refining	
 P30010 Knife Edge Premounted High Shine	 P30001 Wheel Premounted High Shine	


Most Popular

NTI® CeraGlaze®

CeraGlaze® is a remarkable diamond impregnated polishing system used to restore a natural high shine to porcelain, eliminating the need for re-firing. CeraGlaze is bonded with natural caoutchouc. Caoutchouc absorbs heat thus risk of damaging zirconia and other high-performance ceramics. Autoclavable. Handpiece sold individually. Friction grip and RA/ latch sold in packs of 3. HP & RA latch recommended speed: 5000-15,000 rpm. FG polishers at 10,000-12,000 rpm and with light pressure.



First, pre-polish with the green CeraGlaze® which removes scratches from the surface. This polisher can also be used in place of a diamond to make minor adjustments. Next, use the blue refining CeraGlaze® to achieve a lustrous shine and prepare the surface for the final high shine step. After refining the restoration with the blue CeraGlaze®, follow with a yellow reglazing CeraGlaze® polisher to achieve a natural "wet-looking" high shine.



Porcelain

PRE-POLISHING	CERA GLAZE® 1 AXIS / NTI								
	HEAD Ø 1/16mm	250	250	055	145	050	030	060	100
LENGTH (mm)	3.0	2.0	15.5	2.5	10.0	6.0	9.0	4.7	
FG shank - 3 pk						P353-3	P355-3		
RA shank - 3 pk					P332-3	P333-3	P335-3	P336-3	
HP shank - 1pk	P301	P310	P344	P343					opt. 5,000 - 7,000 rpm

Porcelain, Precious / Semi-Precious Metals

REFINING	CERA GLAZE® 2 AXIS / NTI									
	HEAD Ø 1/16mm	250	250	055	150	145	050	030	060	100
LENGTH (mm)	3.0	2.0	15.5	1.0	2.5	10.0	6.0	9.0	4.7	
FG shank - 3 pk							P3053-3	P3055-3		
RA shank - 3 pk							P3032-3	P3033-3	P3035-3	P3036-3
HP shank - 1pk	P3001	P3010	P3044	P3042	P3043					opt. 5,000 - 7,000 rpm

Porcelain, Precious / Semi-Precious Metals

REGLAZE / HIGH SHINE	CERA GLAZE® 3 AXIS / NTI									
	HEAD Ø 1/16mm	250	250	055	150	145	050	030	060	100
LENGTH (mm)	3.0	2.0	15.5	1.0	2.5	10.0	6.0	9.0	4.7	
FG shank - 3 pk							P30053-3	P30055-3		
RA shank - 3 pk							P30032-3	P30033-3	P30035-3	P30036-3
HP shank - 1pk	P30001	P30010	P30044	P30042	P30043					opt. 5,000 - 7,000 rpm

DiaLite is a trademark of Brasseler USA.

Adjusting : Porcelain

HEAD Ø 1/16mm	220	220	170	050	060	045	050	050	
LENGTH (mm)	3.0	Knife Edge	Knife Edge	10	10	16.3	16	12	
Unmounted - 10 pk	P0301-220	P0310-220	P0311-170	P0352	P0355		P0330-045	P0351	P0384
RA shank - 10 pk									
HP shank - 10 pk									

opt. 15,000 rpm

Refining : Porcelain

HEAD Ø 1/16mm	220	220	050	060	045	050	050
LENGTH (mm)	3.0	Knife Edge	10	10	16.3	16	12
Unmounted - 10 pk	P0306-220	P0316-220	P0362	P0365			
RA shank - 10 pk					P0340-045	P0361	P0394
HP shank - 10 pk							

opt. 10,000 rpm

Super High Shine : Porcelain

HEAD Ø 1/16mm	220	220	050	050	150	150	055
LENGTH (mm)	3.0	Knife Edge	10	16	2.0	2.5	15.5
Unmounted - 10 pk	P0321-220	P0322-220					
RA shank - 10 pk			P0372				
HP shank - 10 pk				P0371	P0373	P0375	P0374

opt. 5000 rpm

■ = Bulk Pack (100) add "G" to end of part number

Most Popular

ProGloss™
Half the Time With No Compromise

The KaVo Kerr ProGloss Polisher is a one step polishing system for composite restorations. Achieve a high surface luster faster in one simple step. The multiple use polishing system, easy to handle and highly time saving, is designed for either wet or dry applications, producing excellent results.


- Advance formulations designed for today's composite materials
- Multiple use - One Step Polisher - Reduces chair time
- Flexible design for control and contact with the polishing surface



HEAD Ø 1/16mm	040	070	095
LENGTH (mm)	9.5	9.4	6.2
RA shank - 12 pk	P8001	P8002	P8003
HP shank - 10 pk			


P8000A - set of 4 each opt. 5,000 rpm

Reference Surface




OptiDisc Coarse/Medium
sRa= 0.56 µm


Treated Surface with Single Step Polisher*




ProGloss™
sRa= 0.18 µm




OptraPol
sRa= 0.32 µm



Komet
sRa= 0.19 µm




Jazz
sRa= 0.22 µm



Enhance & PoGo**
sRa= 0.56 µm

*Data on file: ProGloss™ - KaVo Kerr · OptraPol - Vivadent · Komet 1 step Polisher - Gebr. Brasseler · Jazz - SS White · Enhance & PoGo - Dentsply
** two step polisher, according to manufacturer instructions for use.




Mandrels and Shank Converters

Converters are fast and safe for turning your FG instrument into HP shanks. The adapter shank is constructed like a snap-on clamping sleeve. Once inserted, the instrument will remain in the adapter until it has become dull. Made from special stainless steel, the adapter together with the instrument can be sterilized in any normally used solution or in the autoclave.

HEAD Ø 1/16mm	130	050	050	050*	080	140	023	FG to HP	FG to RA
RA shank - 6 pk							RA-SNAP		M024
HP shank - 6 pk	301L	303-050	305-050	305-050LH	305-080	305-14		M022	

*Left Hand Screw Mandrel



UltraGloss™

- Metal free polisher surface prevents damage to restoration
- Flexible design allows polisher to conform to structures
- Super thin polishing disks create natural contours
- Standard color coding for simplified selection



HEAD Ø 1/16mm	120	080	120	080	120	080	120	080
Grit	Coarse	Coarse	Medium	Medium	Fine	Fine	Superfine	Superfine
Surface	Single Sided	Single Sided	Double Sided	Double Sided	Double Sided	Double Sided	Double Sided	Double Sided
Unmounted - 30 pk	DSC-C120	DSC-C080	DSC-M120	DSC-M080	DSC-F120	DSC-F080	DSC-SF120	DSC-SF080

opt. 10,000 - 12,000 rpm



HiLuster PLUS Polishing System

- Maximum flexibility.
- Designed to yield high gloss and superior esthetics.
- Smooth surfaces. Helps prevent plaque buildup and discoloration.
- No polishing paste required. No mess. Improves visibility to work area.
- Durable. Autoclavable. Reusable multiple times.

Gloss PLUS Polisher Refills (pkg of 6) Cup	2653
Gloss PLUS Polisher Refills (pkg of 6) Disc	2654
Gloss PLUS Polisher Refills (pkg of 6) Flame	2651
Gloss PLUS Polisher Refills (pkg of 6) Minipoint	2652
Polishing Kit	210
PLUS Dia Polisher Refills (pkg of 6) Cup	2663
PLUS Dia Polisher Refills (pkg of 6) Disc	2664
PLUS Dia Polisher Refills (pkg of 6) Flame	2661
PLUS Dia Polisher Refills (pkg of 6) Minipoint	2662
HiLuster PLUS Kit	2600



Occlubrush Refills

With Occlubrush, you can finish and polish restorations to a high gloss with one instrument.

- Available in three shapes (regular cup, small cup and point) designed for polishing occlusal and interproximal surfaces.
- Ideal for posterior composites.

Regular Cup (10)	2510
Small Cup (10)	2511
Point (10)	2512

Most Popular



OptiDisc™


The OptiDisc finishing and polishing system gives your restorations a final polish similar to that of natural dentition. With a three-grit system instead of the traditional four, you can achieve a higher surface gloss with fewer steps than with competitive products. OptiDiscs are translucent for enhanced working surface visibility. KaVo Kerr also offers OptiDisc Extra-Coarse for contouring composite restorative materials prior to finishing and polishing.

Features:


- Translucent disc. Enhances work surface visibility.
- Flush-mounted mandrel. Protects neighboring teeth and tissue.
- Identifiable color-coding system. Color indicates grit type and abrasive side.
- Three-step polishing system. Achieves natural high gloss with fewer steps.
- Contouring disc. Extra-Coarse grit available for contouring prior to finishing and polishing.

Assorted Kit	33475
Contains 30 discs each of: 9.6 mm and 12.6 mm – Extra-Fine 9.6 mm and 12.6 mm – Fine, 9.6 mm and 12.6 mm – Coarse/Medium 9.6 mm and 12.6 mm – Extra-Coarse 3 Mandrels	

Precise




Flexible



Transparent



Safe



3M™ SofLex™ to KaVo Kerr OptiDisc Conversions

3M SofLex™

SofLex wXT	SofLex	Description
	1980	Contouring and Polishing Discs Kit, 1980
2382C	1982C	Disc Coarse - Refill, 1/2" Diameter
2382M	1982M	Disc Medium - Refill, 1/2" Diameter
2382F	1982F	Disc Fine - Refill, 1/2" Diameter
2382SF	1982SF	Disc Super Fine - Refill, 1/2" Diameter
2381C	1981C	Disc Coarse - Refill, 3/8" Diameter
2381M	1981M	Disc Medium - Refill, 3/8" Diameter
2381F	1981F	Disc Fine - Refill, 3/8" Diameter
2381SF	1981SF	Disc Super Fine - Refill, 3/8" Diameter

OptiDisc


Order #	Description
33475	Assorted Kit
33483	Refills (pkg of 80) 12.6 mm, Extra-Coarse
33482	Refills (pkg of 80) 12.6 mm, Coarse/Medium
33481	Refills (pkg of 80) 12.6 mm, Fine
33480	Refills (pkg of 80) 12.6 mm, Extra-Fine
33486	Refills (pkg of 80) 15.9 mm, Coarse/Medium
33487	Refills (pkg of 80) 15.9 mm, Extra-Coarse
33484	Refills (pkg of 80) 15.9 mm, Extra-Fine
33485	Refills (pkg of 80) 15.9 mm, Fine
33479	Refills (pkg of 80) 9.6 mm, Extra-Coarse
33478	Refills (pkg of 80) 9.6 mm, Coarse/Medium
33477	Refills (pkg of 80) 9.6 mm, Fine
33476	Refills (pkg of 80) 9.6 mm, Extra-Fine
33488	Mandrels (pkg of 5)

Sof-Lex is a trademarks of 3M Company


DiaGloss Diamond-Impregnated Polisher

The DiaGloss Diamond-Impregnated Polisher is bonded with natural caoutchouc. Caoutchouc absorbs heat thus reducing risk of overheating and deforming composite filings.

Pre-Polishing : Composite

DIAGLOSS®					
	HEAD Ø 1/16mm	050	030	060	100
	LENGTH (mm)	10	6.0	10	1.0
	FG shank - 3 pk		P1953-3		
	RA shank - 3 pk	P1932-3	P1933-3	P1935-3	P1938-3
HP shank - 1 pk					
opt. 5,000 - 7,000 rpm					


High Shine : Composite

DIAGLOSS®					
	HEAD Ø 1/16mm	050	030	060	100
	LENGTH (mm)	10	6.0	10	1.0
	FG shank - 3 pk		P19053-3		
	RA shank - 3 pk	P19032-3	P19033-3	P19035-3	P19038-3
HP shank - 1 pk					
opt. 5,000 - 7,000 rpm					

PDQ2® Composite Polishing



Available in set of 5 each - LS-510

PDQ2™ POLISHERS					
	HEAD Ø 1/16mm	050	030	060	100
	LENGTH (mm)	10	6.0	10	1.0
	FG shank - 10 pk				
	RA shank - 20 pk	P3132	P3133	P3135	P3138
HP shank - 1 pk					
opt. 5,000 - 7,000 rpm					

DiaSheen™

Polishing Paste designed for Composites, Porcelain, Metals & Acrylic or PMMA. DiaSheen is designed to give a high gloss to today's composites, porcelains, ceramics, PMMA and metals. DiaSheen is ideal for use with stiff bristle brushes, felts or other rotary type applications.



- Fine 3 gram tub
- Fine-Soft 3 gram tub
- Medium 3 gram tub

- PP-FINE
- PP-FINESOFT
- PP-MEDIUM

Natural Bristle Brushes

These natural bristle brushes are best for the polishing and fine tuning the post glazed surfaces of occlusal grooves and concave anatomies. Can be used with or without polishing paste.



- Stiff, 0.6875" 12/pack
- Stiff, 0.6875" 144/pack
- Soft, 0.6875" 12/pack
- Soft, 0.6875" 144/pack

- 11-STIFF
- 11-STIFF-G
- 11-SOFT
- 11-SOFT-G

Rag / Muslin Wheel

Ideal for cleaning and polishing of denture acrylic material, cleaning interproximal spaces after denture deflasking, and for general laboratory applications.



- 1 inch / 25.4 mm rag wheel
- RAG-1

PDQ® Diamond-Impregnated Polishing Brush

The PDQ® Diamond-Impregnated Polishing Brush incorporates micro-diamond impregnated bristles for achieving the desired high shine on the new "hard setting" composites. The PDQ® Composite Polishing Brush saves the dental professional valuable chair time. Autoclavable.


Available in set of 5 each - LS-508

RA SHANK				
	SHAPE	Cup	Small Cup	Point
	RA shank - 10 pk	P0800	P0801	P0802
opt. 5,000 rpm max				




☐ = Bulk Pack (100) add "G" to end of part number


Adjusting : Acrylics, Alloys, Precious / Semi-Precious Metals

BLUE / SILICONE						
	HEAD Ø 1/16mm	220	180	220	170	055
	LENGTH (mm)	3.0	4.0	Knife Edge	Knife Edge	16.3
	Unmounted - 10 pk	P0401-220	P0402-180	P0410-220	P0411-170	
	RA shank - 10 pk					P0440-045
HP shank - 10 pk						
opt. 15,000 rpm						


Pre-Polishing : Alloys, Amalgam, Composite, Gold

BROWN / SILICONE										
	HEAD Ø 1/16mm	220	050	060	085	060	050	050	150	145
	LENGTH (mm)	3.0	10	6.0	7.0	10	16	12	2.0	2.5
	Unmounted - 10 pk	P0001-220								
	FG shank - 10 pk			P0053-031		P0055-065				
RA shank - 10 pk		P0032-047	P0033-031	P0034-045	P0035-065					
HP shank - 10 pk			P0046			P0041	P0044	P0042-150	P0043-150	
opt. 5,000 - 7,000 rpm										

High Shine : Alloys, Amalgam, Composite, Gold

GREEN / SILICONE										
	HEAD Ø 1/16mm	220	050	060	085	060	050	050	150	145
	LENGTH (mm)	3.0	10	6.0	7.0	10	16	12	2.0	2.5
	Unmounted - 10 pk									
	FG shank - 10 pk			P0153-031		P0155-065				
RA shank - 10 pk		P0132-047	P0133-031	P0134-045	P0135-065					
HP shank - 10 pk	P0101-220		P0146			P0141	P0144	P0142-150	P0143-150	
opt. 10,000 - 15,000 rpm										

Gold, Precious / Semi-Precious Metals

GOLDSTAR / SILICONE										
	HEAD Ø 1/16mm	220	055	050	220	055	050	220	055	050
	LENGTH (mm)	4.0	15.5	12	4.0	15.5	12	4.0	15.5	12
	Unmounted - 10 pk	P1801			P1811			P1812	P1821	
	RA shank - 10 pk			P18032				P18132		P18232
HP shank - 10 pk		P18044			P18144			P18244		
opt. 20,000 rpm										
opt. 15,000 rpm										
opt. 10,000 rpm										

Assortment pack (5 of each) of 22mm unmounted wheels - P1899

Shammy Bright™ Polishing System

Polishing Disc and Paste for Composites and Porcelain

A complete set of pastes and polishing discs for the final polishing of composite and porcelain restorations.

- 2g syringe, 3-micron diamond paste
- 5g syringe, Sub-micron a luminous oxide paste
- 0.5 inch, Micro-cloth polishing disc (25 per pack)
- 0.3125 inch, Micro-cloth polishing disc (25 per pack)

- PPF-COM
- PPF-POR
- DSC-LD
- DSC-SD

