



Midwest® Carbide Burs



Round

Useful for opening and shaping cavity preparations. Also used for root canal access and caries removal. Smaller sizes facilitate conservative preparations.

Round

Shape		1/4	1/2	1	2	3	4	5	6	7	8	10
Diameter mm		0.5	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.1	2.3	2.7
Length mm		0.5	0.6	0.8	1.0	1.2	1.4	1.6	1.8	2.1	2.3	2.7
FG	100 Pk Bag	386201	386202	386203	386204	386268	386205	386269	386206	386274	386207	
	100 Pk Box	385201	385202	385203	385204	385268	385205	385269	385206	385274	385207	
	10 Pk	389201	389202	389203	389204	389268	389205	389269	389206	389274	389207	
	Sterile 25 Pk	314201	314202	314203	314204		314205		314206		314207	
FGSS	100 Pk Bag	386401	386402	386403	386404		386405		386406			
	100 Pk Box	385401	385402	385403	385404		385405					
	10 Pk	389401	389402	389403	389404		389405					
	Sterile 25 Pk				314404							
LA	100 Pk Bag	386101	386102	386103	386104		386105	386169	386106		386107	
	100 Pk Box	385101	385102	385103	385104		385105	385169	385106		385107	
	10 Pk	389101	389102	389103	389104		389105	389169	389106		389107	
	Sterile 25 Pk		314102	314103	314104		314105		314106		314107	
HP	100 Pk Bag	386301	386302	386303	386304		386305		386306		386307	
	10 Pk	389301	389302	389303	389304		389305		389306		389307	389336
	Sterile 10 Pk				319304		319305					



Inverted Cone

Ideal for creating mechanical retentions and undercuts when contouring cavity preparations for amalgam restorations.

Inverted Cone

Shape		33 1/2	34	35	36	37	39
Diameter mm		0.6	0.9	1.0	1.0	1.2	1.2
Length mm		0.4	0.6	0.7	0.8	0.8	1.3
FG	100 Pk Bag	386209	386210	386211	386212	386213	386215
	100 Pk Box	385209	385210	385211	385212	385213	385215
	10 Pk	389209	389210	389211	389212	389213	389215
	Sterile 25 Pk	314209	314210	314211	314212	314213	
FGSS	100 Pk Bag	386409	386410	386411		386413	
	100 Pk Box	385409	385410	385411		385413	
	10 Pk	389409	389410	389411		389413	
	Sterile 25 Pk			314411			
LA	100 Pk Bag	386109	386110	386111	386112	386113	
	100 Pk Box	385109	385110	385111	385112	385113	
	10 Pk	389109	389110	389111	389112	389113	
	Sterile 25 Pk		314110	314111			
HP	100 Pk Bag	386309				386313	386315
	10 Pk	389309				389313	389315




Midwest® Operative Carbide Burs



Pear

Designed for easy initial entry and a smooth contoured shape that eliminates sharp corners. Can also be used to create convergent lateral walls or mechanical retention.

Pear




	Shape	329	330	331	331-L	332	332-L	333-L
	Diameter mm	0.6	0.8	1.0	1.0	1.2	1.2	1.4
	Length mm	1.6	1.6	1.6	4.1	1.6	4.1	4.1
FG	100 Pk Bag	386369	386261	386262	386264	386263	386265	386266
	100 Pk Box	385369	385261	385262	385264	385263	385265	385266
	10 Pk	389369	389261	389262	389264	389263	389265	389266
	Sterile 25 Pk	314269	314261	314262	314264	314263		
FGSS	100 Pk Bag	386408	386461	386462	386464	386463		
	100 Pk Box	385408	385461	385462	385464	385463		
	10 Pk	389408	389461	389462	389464	389463		
	Sterile 25 Pk		314461	314462				
LA	100 Pk Bag		386161					
	10 Pk		389161	389162				



Pear Fine Crosscut


	Shape		1930	1931		1932		
FG	100 Pk Bag			386559				
	100 Pk Box			385559				
	10 Pk		389551	389559		389552		
	2 Pk		388551	388559		388552		

Long Pear



	Shape	245	246
	Diameter mm	0.9	1.2
	Length mm	2.7	3.7
FG	100 Pk Bag	386371	386372
	100 Pk Box	385371	385372
	10 Pk	389371	389372
	Sterile 25 Pk	314271	
FGSS	100 Pk Bag	386471	
	100 Pk Box	385471	
	10 Pk	389471	
	Sterile 25 Pk	314471	

All Purpose




	Shape	256	257
	Diameter mm	0.9	1.0
	Length mm	2.7	3.7
FG	100 Pk Bag	386287	386288
	100 Pk Box	385287	385288
	10 Pk	389287	389228
	Sterile 25 Pk		
FGSS	100 Pk Bag		
	100 Pk Box		
	10 Pk		
	Sterile 25 Pk		

Inlay Prep.



	Shape	271
	Diameter mm	1.2
	Length mm	2.7
FG	100 Pk Bag	386373
	100 Pk Box	385373
	10 Pk	389373
	Sterile 25 Pk	
FGSS	100 Pk Bag	
	100 Pk Box	
	10 Pk	
	Sterile 25 Pk	

End Cutting




	Shape	956	957
	Diameter mm	0.9	1.0
	Length mm	3.2	3.7
FG	100 Pk Bag	386276	386278
	100 Pk Box	385276	
	10 Pk		389278
	Sterile 25 Pk		
FGSS	100 Pk Bag		
	100 Pk Box		
	10 Pk		
	Sterile 25 Pk		

Straight Flat End

Ideal for creating flat or parallel surfaces; can be used to create flat floors with square shoulders.


Extended length heads and crosscut features also available.



Straight Dome End

Designed for easier entry than a flat end bur. Creates a straight side surface with smooth rounded corners. Extended length heads and crosscut versions also available.

STRAIGHT FLAT END



Diameter mm	0.9	0.9	1.0	1.0	1.2	1.2	1.6	1.6
Length mm	3.2	4.9	3.7	4.9	3.7	4.9	4.2	5.7

Straight Flat End



	Shape	56	57	57-L	58		60-L
FG	100 Pk Bag	386216	386217	386240	386218		
	100 Pk Box	385216	385217	385240	385218		
	10 Pk	389216	389217	389240	389218		
	Sterile 25 Pk	314216	314217		314218		
FGSS	100 Pk Bag	386416	386417				
	100 Pk Box	385416	385417				
	10 Pk	389416	389417				
LA	100 Pk Bag	386116					
	100 Pk Box	385116					
HP	100 Pk Bag						386320
	10 Pk						

Straight Flat End Crosscut




	Shape	556	556-L	557	557-L	558	558-L	560
FG	100 Pk Bag	386226	386290	386227	386291	386228	386292	
	100 Pk Box	385226	385290	385227	385291	385228		
	10 Pk	389226	389290	389227	389291	389228	389292	
	Sterile 25 Pk	314226		314227	314251	314228		
FGSS	100 Pk Bag	386426		386427	386491	386428		
	100 Pk Box	385426		385427	385491	385428		
	10 Pk	389426		389427	389491	389428		
	Sterile 25 Pk	314426		314427				
LA	100 Pk Bag	386126		386127		386128		
	100 Pk Box	385126		385127				
	10 Pk	389126		389127		389128		
HP	100 Pk Bag			386327		386328		
	10 Pk			389327		389328		389330
	Sterile 10 Pk			319327				

Straight Flat End Fine Crosscut



	Shape		2057				
FG	100 Pk Bag		386560				
	100 Pk Box		385560				
	10 Pk		389560				
	2 Pk		388560				

STRAIGHT DOME END



Diameter mm	0.9	1.0	1.2	1.2
Length mm	3.2	3.7	3.7	4.9

Straight Dome End



	1156	1157	1158	
	386255	386256	386257	
	385255	385256	385257	
	389255	389256	389257	
	314255	314256	314257	
	386455	386456	386457	
	385455	385456	385457	
	389455	389456	389457	
		386356		
		389356		

Straight Dome End Crosscut




	1556	1557	1558	
	386258	386259	386260	
	385258	385259	385260	
	389258	389259	389260	
	314258	314259	314260	
	386458	386459	386460	
	385458	385459	385460	
	389458	389459	389460	
		386359		
		389359	389360	

Straight Dome End Fine Crosscut



	1956	1957	1958	1958-L
			386558	
			385558	
	389556	389557	389558	389555
	388556	388557	388558	388555

TAPER FLAT END



Taper Flat End

	Shape	167	168	169	169-L	170	170-L	171	171-L			
FG	100 Pk Bag		386370	386221	386245	386222	386246	386223	386247			
	100 Pk Box		385370	385221	385245	385222	385246	385223	385247			
	10 Pk	389801		389221	389245	389222	389246	389223	389247			
	Sterile 25 Pk			314221	314245	314222	314246					
FGSS	100 Pk Bag			386419		386422	386446					
	100 Pk Box			385419		385422	385446					
	10 Pk			389419		389422	389446					
LA	100 Pk Bag				386145							
	100 Pk Box				385145							
	10 Pk			389121	389145	389122						



Taper Flat End Crosscut

	Shape			699	699-L	700	700-L	701	701-L	702	703	703-L
FG	100 Pk Bag			386231	386248	386232	386249	386233	386250	386234	386235	
	100 Pk Box			385231	385248	385232	385249	385233	385250	385234	385235	
	10 Pk			389231	389248	389232	389249	389233	389250	389234	389235	
	Sterile 25 Pk			314231		314232		314233	314250	314234		
FGSS	100 Pk Bag			386431		386432		386433				
	100 Pk Box			385431		385432		385433				
	10 Pk			389431		389432		389433				
LA	100 Pk Bag					386132		386133				
	100 Pk Box					385132		385133				
	10 Pk					389132		389133		389134		
HP	100 Pk Bag					386332		386333	386350	386334	386335	386367
	10 Pk					389332		389333	389350	389334	389335	389367
	Sterile 10 Pk							319350				



Taper Flat End Fine Crosscut

	Shape			2069		2070						
FG	10 Pk			389568								
	2 Pk			388568		388569						





Taper Dome End

Can be used to create divergent sides with rounded corners. Designed for easy initial entry with the rounded end. Use crosscut versions for more efficient gross cutting.

TAPER DOME END

Diameter mm	0.9	0.9	1.0	1.0	1.2	1.2	1.6	1.6	2.1	2.1
Length mm	3.2	4.9	3.7	4.9	3.7	4.9	4.1	4.9	4.5	5.7



Taper Dome End

	Shape	1169	1169-L	1170	1170-L	1171	1171-L	1172			
FG	100 Pk Bag	386271	386242	386272	386243	386273	386244	386294			
	100 Pk Box		385242	385272	385243	385273	385244	385294			
	10 Pk	389271	389242	389272	389243	389273	389244	389294			
	Sterile 25 Pk					314273					
LA	100 Pk Bag					386189		386194			
	100 Pk Box					385189		385194			
	10 Pk					389189		389194			
HP	100 Pk Bag					386346		386394			
	10 Pk					389346		389394			



Taper Dome End Crosscut

	Shape			1700		1701		1702	1702-L	1703	1703-L
FG	100 Pk Bag			386282		386284		386286			
	100 Pk Box			385282		385284		385286			
	10 Pk			389282		389284		389286			
HP	10 Pk							389386	389380	389385	389376




Taper Dome End Fine Crosscut

	Shape			1970		1971					
FG	10 Pk			389561		389572					
	2 Pk			388561		388572					

SPECIALTY METAL CUTTING BURS


Exceptional metal cutting performance from our unique blade geometry, designed for optimum endurance and efficiency.

Straight Dome End CrissCross Cut





	Shape	2157	2158	2170	2171
	Diameter mm	1.0	1.2	1.0	1.2
	Length mm	3.7	3.7	3.7	3.7
FG	100 Pk Bag		386562		
	100 Pk Box		385562		
	10 Pk	389563	389562	389564	389565
	2 Pk	388563	388562	388564	




Taper Dome End CrissCross Cut




Midwest® Oral Surgical Carbide Burs

Superior, precision ground cutting edges for efficient and consistent cutting performance. Ideal for cutting through bone and sectioning teeth.

		Round								Inverted Cone			
													
													
		Shape	1/4	1/2	1	2	4	6	8	10	33 1/2	37	39
		Diameter mm	0.5	0.6	0.8	1.0	1.4	1.8	2.3	2.7	0.6	1.4	1.8
		Length mm	0.5	0.6	0.8	1.0	1.4	1.8	2.3	2.7	0.4	0.8	1.3
FGOS	100 Pk Box						385605	385606	385607				
	10 Pk			389603									
	2 Pk				388604	388605	388606	388607					
	Sterile 10 Pk		314602	314603	314604	314605	314606	314607					
LAOS	100 Pk Box						385905	385906	385907				
	10 Pk			389992									
	2 Pk				388904	388905	388906	388907					
	Sterile 10 Pk				314904	314905	314906	314907					
HP	100 Pk Bag	386301	386302	386303	386304	386305	386306	386307			386309	386313	386315
	10 Pk	389301	389302	389303	389304	389305	389306	389307	389336		389309	389313	389315


		Straight Flat End		Straight Flat End Crosscut					Straight Dome End		Straight Dome End Crosscut	
												
												
		Shape	57	60-L	556	557	558	559	560	1157	1557	1558
		Diameter mm	1.0	1.6	0.9	1.0	1.2	1.4	1.6	1.0	1.0	1.2
		Length mm	3.7	5.7	3.2	3.7	3.7	4.2	4.2	3.7	3.7	3.7
FGOS	10 Pk				389626	389627						
	2 Pk	388617			388626	388627	388628	388629			388659	388660
	Sterile 10 Pk					314627	314628	314629				
LAOS	10 Pk					389927	389928					
	2 Pk					388927						388960
	Sterile 10 Pk					314927						
HP	100 Pk Bag		386320			386327	386328			386356	386359	386360
	10 Pk					389327	389328		389330	389356	389359	389360

Taper Flat End Crosscut




	Shape	700	701	701-L	702	702-L	703	703-L
	Diameter mm	1.0	1.2	1.2	1.6	1.6	2.1	2.1
	Length mm	3.7	3.7	4.9	4.1	4.9	4.5	5.7
FGOS	100 Pk Bag					386651		
	100 Pk Box					385651		
	10 Pk	389632				389651		
	2 Pk	388632	388633		388634		388635	
	Sterile 10 Pk	314632	314633		314634		314635	
LAOS	10 Pk	389932			389934			
	2 Pk		389933				389935	
HP	100 Pk Bag	386332	386333	386350	386334		386335	386367
	10 Pk	389332	389333	389350	389334		389335	389367
HPOS	10 Pk				389609		389735	


Taper Dome End




Taper Dome End Crosscut



Taper Dome End Long



Taper Safe End



	Shape	1171	1172	1702	1702-L	1703	1703-L	151-Z	151L-Z	152
	Diameter mm	1.2	1.6	1.6	1.6	2.1	2.1	1.6	1.6	1.6
	Length mm	3.7	4.1	4.1	5.3	4.5	5.7	11.0	11.0	9.0
FG	Sterile 10 Pk							319251		319224
FGOS	10 Pk					389685				
	2 Pk			388686						
	Sterile 10 Pk							319252		
LAOS	10 Pk			389986		389985				
HP	100 Pk Bag	386346	386394							
	10 Pk	389346	389394	389386	389380	389385	389376			
HPOS	10 Pk					389785				

Lindemann




Lindemann Crosscut



	Shape	161	162
	Diameter mm	1.6	1.6
	Length mm	10.0	10.0
FGOS	Sterile 10 Pk	314613	314614

Midwest® Trimming & Finishing Carbide Burs




Round

May be used for contouring and finishing of restorations, root planning or contouring and smoothing of bone.



Egg

Geometrically suited to cutting concave surfaces: occlusals and anterior lingual.



Needle

Narrow tapered end for placing secondary anatomy on restorations.



Interproximal


An ideal curvature for interproximal trimming and fine finishing on anterior and posterior restorations.



Taper Flat End

Tapered burs with flat ends, ideal for finishing composite restorations, available in 12 and 30 blade options.

Round




12 Blade	Shape	7002	7003	7004	7006	7008	7009
	Diameter mm	1.0	1.2	1.4	1.6	2.3	2.5
	Length mm	1.0	1.2	1.4	1.6	2.3	2.5
FG	100 Pk Bag						
	100 Pk Box						
	10 Pk	389501	389502	389503	389504	389505	389506
	2 Pk	388501	388502	388503	388504	388505	388506
	Sterile 10 Pk	319504					



30 Blade	Shape	9004	9006	9008
FG	2 Pk	388591	388592	388593
	Sterile 10 Pk			

Egg




12 Blade	Shape	7404	7406	7408
	Diameter mm	1.4	1.8	2.3
	Length mm	3.0	3.4	3.6
FG	100 Pk Bag	386522	386523	386524
	100 Pk Box	385522	385523	385524
	10 Pk	389522	389523	389524
	2 Pk	388522	388523	388524
	Sterile 10 Pk	319522	319523	319524



30 Blade	Shape	9406
FG	2 Pk	388584
	Sterile 10 Pk	319223

Bullet




12 Blade	Shape	7801	7802	7803
	Diameter mm	0.9	1.0	1.2
	Length mm	2.8	2.8	3.0
FG	100 Pk Bag	386507	386508	386509
	100 Pk Box		385508	385509
	10 Pk	389507	389508	389509
	2 Pk	388507	388508	388509
	Sterile 10 Pk		319508	319509



30 Blade	Shape	9803
FG	2 Pk	388599
	Sterile 10 Pk	

Needle



12 Blade	Shape	7901	7902	7903	
	Diameter mm	0.9	1.0	1.2	1.4
	Length mm	3.2	3.6	3.8	5.2
FG	100 Pk Bag	386510	386511	386512	
	100 Pk Box	385510	385511	385512	
	10 Pk	389510	389511	389512	
	2 Pk	388510	388511	388512	
	Sterile 10 Pk	319510	319511	319512	



30 Blade	Shape	9903	9904
FG	2 Pk	388597	388598
	Sterile 10 Pk	319597	

Interproximal



12 Blade	Shape	7103
	Diameter mm	1.2
	Length mm	3.4
FG	100 Pk Bag	
	100 Pk Box	
	10 Pk	389546
	2 Pk	388546
	Sterile 10 Pk	



30 Blade	Shape	9103
FG	2 Pk	388581
	Sterile 10 Pk	319546

Flame




12 Blade	Shape	274-016	7104	7106	7108
	Diameter mm	1.6	1.4	1.8	2.3
	Length mm	3.5	3.2	3.6	4.4
FG	10 Pk		389513	389514	389515
	2 Pk		388513	388514	388515
	Sterile 10 Pk	319576	319513	319514	

Flame Long



12 Blade	Shape	48L-010	48L-012
	Diameter mm	1.0	1.2
	Length mm	8.0	8.0
FG	10 Pk		
	2 Pk		
	Sterile 10 Pk	319547	319548

Inverted Taper



12 Blade	Shape	7303	7304
	Diameter mm	1.2	1.4
	Length mm	4.1	4.1
FG	10 Pk	389519	389520
	2 Pk	388519	388520

		Taper Flat End		Taper Cone			Straight Dome End		Straight Flat End		Chamfering	
												
												
12 Blade	Shape	7702	7713	7204	7205	7206			7572		7345	
	Diameter mm	1.0	1.2	1.4	1.6	1.8	0.9		1.0		1.5	
	Length mm	4.9	4.9	8.8	8.8	8.8	3.4		5.0		—	
FG	10 Pk	389528	389529	389516	389517	389518			389525		389545	
	2 Pk	388528	388529	388516	388517	388518			388525		388545	
												
30 Blade	Shape		9713				9561		9572			
FG	2 Pk		388596				388594		388595			

		Taper T-Series			Taper Long & Extra Long		
							
							
12 Blade	Shape	7114 T-1	7214 T-2	7714 T-3	7642	7664	7675
	Diameter mm	0.9	1.2	1.4	1.0	1.5	1.6
	Length mm	2.9	5.1	7.6	5.7	8.6	8.0
FG	100 Pk Bag					386537	
	100 Pk Box					385537	
	10 Pk	389901	389902	389903	389536	389537	
	2 Pk	388901	388902	388903	388536	388537	
					319568	319569	319570
							
30 Blade	Shape		9214 T-2	9714 T-3	9642		
FG	2 Pk		388583	388582	388585		
	Sterile 10 Pk		319501				

		OF Series • Safe End		
				
				
12 Blade	Shape	OF1	OF2	OF3
	Diameter mm	2.3	0.9	0.7
	Length mm	4.3	3.2	3.2
FG	Sterile 10 Pk	319588	319589	319590

OF1 Football Safe End
OF2, OF3 Taper Safe End

		EF Series • Taper Safe End			
					
					
10 Blade	Shape	EF3	EF4	EF6	EF9
	Diameter mm	0.8	1.0	1.4	1.4
	Length mm	3.0	4.0	6.0	9.0
FG	Sterile 10 Pk	319579	319580	319581	319582
30 Blade	Shape				EF9-UF
FG	Sterile 10 Pk				319583

SPECIALTY

Crown Preparation Series Metal Cutting • Fine Crosscut (FCC)

							
							
12 Blade	Shape	1931	1958	379-023	847-016	856-016	856-018
	Diameter mm	1.0	1.2	2.3	1.6	1.6	1.8
	Length mm	1.9	4.0	4.3	8.0	8.0	8.0
FG	Sterile 25 Pk	314277	314278	314279	314247	314286	314288

1931 Pear FCC
1958 Straight Dome End FCC
379 Football FCC
847 Taper Flat End Mod. FCC
856 Taper Dome End FCC


Midwest® Esthetic Finishing Carbide Burs

- Multiple precision ground edges for exceptionally smooth surface finish.
- 21 anatomically adapted shapes to meet every finishing challenge.
- Offered in regular and short shanks with available safe end tips.

			Egg			Bullet			
									
									
									
			Shape	7404B	7406B	7408B	7801	7802	7803
			Diameter mm	1.4	1.8	2.3	0.9	1.0	1.2
			Length mm	3.0	3.4	3.6	2.8	2.8	3.0
FG	100 Pk			376117		376119	376132	376133	376134
	10 Pk			379117	379118	379119	379132	379133	379134
	2 Pk			378117	378118	378119	378132	378133	378134
FGSS	10 Pk			379114	379115	379116	379136	379137	379138
	2 Pk			378114	378115	378116	378136	378137	378138
									
			Shape	9404B	9406B	9408B			9803
FG	2 Pk			378123	378124	378125			378135
FGSS	2 Pk			378120	378121	378122			378139


			Taper			Taper Safe End			Interproximal	
										
										
			Shape	7611	7612	7613	7114S	7214S	7714S	7103
			Diameter mm	1.0	1.2	1.4	0.9	1.2	1.4	1.2
			Length mm	3.9	5.9	8.8	2.9	5.1	7.6	3.4
FG	10 Pk			379102	379107	379111	379100	379105	379109	379130
	2 Pk			378102	378107	378111	378100	378105	378109	378130
FGSS	10 Pk			379128	379140	379129	379103	379126	379112	
	2 Pk			378128	378140	378129		378126	378112	
										
			Shape				9114S	9214S	9714S	9103
FG	2 Pk						378101	378106	378110	
FGSS	2 Pk						378104	378127	378113	378131





Round

Equally well suited to rapid caries removal, restoration removal, endo access and endo access through a crown. MultiPrep® round burs are truly versatile, multi-purpose burs.



Pear


MultiPrep® Pear shape burs offer excellent multipurpose cutting capabilities for restoration removal, cavity preparation or endo access. Available in Pear and Long Pear shapes.



Straight Flat End


MultiPrep® burs are an ideal choice for old restoration removal and other high performance cutting duties. Use MultiPrep® burs as an upgrade to the standard burs (i.e. 557).

Round Crosscut (CC)




Shape		2	4	6	8
Diameter mm		1.0	1.4	1.8	2.3
Length mm		1.0	1.4	1.8	2.3
FG	100 Pk Bag	396204	396205	396206	396207
	100 Pk Box	395204	395205	395206	395207
	10 Pk	399204	399205	399206	399207
FGSS	100 Pk Bag				
	100 Pk Box				
	10 Pk				
FGOS	100 Pk Bag			396606	
	100 Pk Box			395606	
	10 Pk			399606	
	2 Pk			398606	

Pear CC



Shape		329	330	331	332
Diameter mm		0.6	0.8	1.0	1.2
Length mm		1.6	1.6	1.6	1.6
FG	100 Pk Bag		396261	396262	396263
	100 Pk Box	395369	395261	395262	395263
	10 Pk	399369	399261	399262	399263
FGSS	100 Pk Bag		396461	396462	
	100 Pk Box		395461	395462	
	10 Pk		399461	399462	
FGOS	100 Pk Bag				
	100 Pk Box				
	10 Pk				
	2 Pk				

Long Pear CC




Shape		245	246
Diameter mm		0.9	1.2
Length mm		2.7	3.7
FG	100 Pk Bag	396371	396372
	100 Pk Box	395371	395372
	10 Pk	399371	399372
FGSS	100 Pk Bag		
	100 Pk Box		
	10 Pk		
FGOS	100 Pk Bag		
	100 Pk Box		
	10 Pk		
	2 Pk		

Straight Flat End CC




Shape		556	557	558
Diameter mm		0.9	1.0	1.2
Length mm		3.2	3.7	3.7
FG	100 Pk Bag	396226	396227	396228
	100 Pk Box		395227	395228
	10 Pk	399226	399227	399228
FGSS	100 Pk Bag		396427	
	100 Pk Box		395427	
	10 Pk		399427	
FGOS	100 Pk Bag		396627	
	100 Pk Box		395627	
	10 Pk		399627	
	2 Pk		398627	

Straight Dome End CC




Shape		1556	1557	1558	1558-L
Diameter mm		0.9	1.0	1.2	1.2
Length mm		3.2	3.7	3.7	4.9
FG	100 Pk Bag	396258	396259	396260	396555
	100 Pk Box	395258	395259	395260	395555
	10 Pk	399258	399259	399260	399555
FGSS	100 Pk Bag			396460	
	100 Pk Box			395460	
	10 Pk			399460	
FGOS	100 Pk Bag				
	100 Pk Box				
	10 Pk				
	2 Pk				

Taper Flat End CC




Shape		699	700	701	702
Diameter mm		0.9	1.0	1.2	1.6
Length mm		3.2	3.7	3.7	4.1
FG	100 Pk Bag	396231		396233	396234
	100 Pk Box	395231	395232	395233	395234
	10 Pk	399231	399232	399233	399234

Taper Dome End CC



Shape		1700	1701
Diameter mm		1.0	1.2
Length mm		3.7	3.7
FG	100 Pk Bag	396282	396284
	100 Pk Box	395282	
	10 Pk	399282	399284

Inverted Cone CC



Shape		34	35	37
Diameter mm		0.9	1.0	1.2
Length mm		0.6	0.7	0.8
FG	100 Pk Bag	396210	396211	396213
	100 Pk Box		395211	395213
	10 Pk	399210	399211	399213