


Safco FG Diamonds

Made from carefully selected natural diamond crystals.
Precision-crafted in Switzerland for long-lasting cutting efficiency.




	Pointed football		Football			Interproximal	Ball
Bur No. / ISO#	368-018	368-023	379-018	379-023	390-014	392-016	801-009
Diameter mm	1.8	2.3	1.8	2.3	1.6	1.6	0.9
Length mm	3.5	5.0	3.5	4.5	3.0	5.0	
X-Coarse ●		#102-2549	#102-1902	#102-2554			
Coarse ●	#102-1876	#102-1905	#102-1914	#102-1928	#102-1931		
Regular	#102-1873	#102-1901	#102-1912	#102-1923	#102-1933	#102-1004	#102-1002
Fine (50µm) ●	#102-1879	#102-1909	#102-1917	#102-1926			
Fine (40µm) ●		#102-1910	#102-1918	#102-1929	#102-1935	#102-1007	
Regular SS		#102-1896					#102-1006

	Ball (con't)					Inverted cone	
Bur No. / ISO#	801-012	801-014	801-016	801-018	801-023	805-012	805-016
Diameter mm	1.2	1.8	2.3	1.8	1.6	1.2	1.6
Length mm		3.5	5.0	3.5	4.5	1.0	1.5
X-Coarse ●			#102-2515	#102-2512	#102-2568		#102-2521
Coarse ●	#102-1017	#102-1028	#102-1035	#102-1046	#102-1057	#102-1064	#102-1075
Regular	#102-1013	#102-1024	#102-1031	#102-1042	#102-1053	#102-1061	#102-1072
Fine (50µm) ●				#102-1050	#102-1059		
Fine (40µm) ●					#102-1062		
Regular SS	#102-1020					#102-1068	

	Inverted cone w/collar	Inverted cone long		Barrel		Amalgam remover	
Bur No. / ISO#	806-012	807-012	807-014	811-033	811L-037	813-016	813-020
Diameter mm	1.2	1.2	1.4	3.3	3.7	1.6	2.0
Length mm	2.5	4.0	5.0	4.5	6.0	1.5	2.5
X-Coarse ●				#102-1116	#102-2545		
Coarse ●	#102-1089	#102-1110	#102-1114	#102-1118	#102-1121		
Regular	#102-1084	#102-1106	#102-1112		#102-1125	#102-1134	#102-1142
Regular SS	#102-1086						

Safco FG Diamonds

Made from carefully selected natural diamond crystals.
Precision-crafted in Switzerland for long-lasting cutting efficiency.



Bur No. / ISO#	Pear				Cylinder flat end		
Bur No. / ISO#	830L-014	830R-010	830R-012	830R-018	835-009	835-010	835-010
Diameter mm	1.4	1.0	1.2	1.8	0.9	1.0	1.9
Length mm	4.5	2.8	3.0	3.0	4.0	3.0	4.0
X-Coarse ●				#102-2526			
Coarse ●			#102-1157	#102-1165	#102-1203	#102-1186	#102-1218
Regular	#102-1144	#102-1149	#102-1154			#102-1181	#102-1215
Fine (50µm) ●				#102-1169		#102-1190	
Regular SS			#102-1160		#102-1209	#102-1192	

Bur No. / ISO#	Cylinder flat end (cont.)					Cylinder round end
Bur No. / ISO#	835-012	835-014	836-012	836-014	837-012	838-009
Diameter mm	1.2	1.4	1.2	1.4	1.2	0.9
Length mm	4.0	4.0	6.0	6.0	8.0	4.0
X-Coarse ●				#102-2508		
Coarse ●	#102-1226	#102-1244	#102-1259	#102-1265	#102-1274	
Regular	#102-1222	#102-1242	#102-1254	#102-1263	#102-1272	#102-1296
Fine (50µm) ●	#102-1231	#102-1246				
Regular SS	#102-1237	#102-1248		#102-1267		
Coarse SS ●				#102-1271		

Bur No. / ISO#	Flat end taper							
Bur No. / ISO#	845-012	845-014	846-018	847-012	847-016	847-018	847KR-016	848-016
Diameter mm	1.2	1.4	1.8	1.2	1.6	1.8	1.6	1.6
Length mm	4.0	4.0	6.0	8.0	8.0	8.0	8.0	10.0
X-Coarse ●			#102-2504					#102-2571
Coarse ●	#102-1304	#102-1316	#102-1327	#102-1336	#102-1365	#102-1385	#102-1410	#102-1416
Regular	#102-1301	#102-1313	#102-1322	#102-1333	#102-1361	#102-1382	#102-1403	#102-1414
Fine (50µm) ●		#102-1319			#102-1367	#102-1388		#102-1420
Regular SS	#102-1310			#102-1338	#102-1376			
Coarse SS ●				#102-1340	#102-1379			


	Round end taper			Tree		Round end taper	
Bur No. / ISO#	850-012	850-016	850-018	852-014	852-016	855-012	855-014
Diameter mm	1.2	1.6	1.8	1.4	1.6	1.2	1.4
Length mm	10.0	10.0	10.0	6.0	8.0	7.0	7.0
X-Coarse ●		#102-2579	#102-2537				
Coarse ●	#102-1430	#102-1438	#102-1446	#102-1459	#102-1465	#102-1484	#102-1509
Regular	#102-1427	#102-1435	#102-1443	#102-1451	#102-1462	#102-1481	#102-1502
Fine (50µm) ●			#102-1449		#102-1468	#102-1486	
Regular SS						#102-1498	
Coarse SS ●						#102-1500	

	Round end taper (con't)						
Bur No. / ISO#	855-016	855-025	856-012	856-014	856-016	856-018	856-025
Diameter mm	1.6	2.5	1.2	1.4	1.6	1.8	2.5
Length mm	6.0	7.0	8.0	8.0	8.0	9.0	8.0
X-Coarse ●		#102-2590	#102-1513		#102-1534	#102-2533	#102-1540
Coarse ●	#102-1476	#102-1521	#102-1515	#102-1526	#102-1525	#102-1537	#102-1543
Regular	#102-1473	#102-1517			#102-1522	#102-1533	
Fine (50µm) ●	#102-1482				#102-1518	#102-1538	
Fine (40µm) ●				#102-1531			#102-1545
Regular SS					#102-1530		
Coarse SS ●					#102-1529		

	Needle						
Bur No. / ISO#	858-010	858-012	858-014	859-010	859-012	859-014	859-016
Diameter mm	1.0	1.2	1.4	1.0	1.2	1.4	1.6
Length mm	7.0	7.0	8.0	10.0	11.0	10.0	10.0
Coarse ●				#102-1567	#102-1573	#102-1586	#102-1597
Regular			#102-1562		#102-1571	#102-1589	#102-1594
Fine (50µm) ●	#102-1550	#102-1553	#102-1566				#102-1600
Fine (40µm) ●			#102-1564	#102-1569	#102-1572	#102-1591	#102-1601
X-fine (15µm) ●			#102-1568		#102-1575		#102-1603
Regular SS		#102-1558					
Coarse SS ●		#102-1560					

Safco FG Diamonds

Made from carefully selected natural diamond crystals.
Precision-crafted in Switzerland for long-lasting cutting efficiency.



	Flame						
Bur No. / ISO#	860-014	861-012	861-014	861-018	862-010	862-014	862-016
Diameter mm	1.4	1.2	1.4	1.8	1.0	1.4	1.6
Length mm	5.0	7.0	8.0	6.5	8.0	8.0	8.0
Coarse ●	#102-1616	#102-1624	#102-1647	#102-1662	#102-1663	#102-1665	#102-1678
Regular	#102-1613	#102-1622	#102-1644	#102-1655		#102-1661	#102-1672
Fine (50µm) ●	#102-1619	#102-1626					
Fine (40µm) ●		#102-1628			#102-1666		
X-fine (15µm) ●		#102-1630					
Regular SS		#102-1637					
Coarse SS ●		#102-1640					

	Flame (con't)			Beveled cylinder	Curettage		
Bur No. / ISO#	863-012	863-014	863-016	877-012	877K-014	878K-012	878K-014
Diameter mm	1.2	1.4	1.6	1.2	1.4	1.2	1.4
Length mm	10.0	10.0	10.0	6.0	6.0	8.0	8.0
X-Coarse ●				#102-1739	#102-1748		
Coarse ●	#102-1687	#102-1694	#102-1701	#102-1735	#102-1744	#102-1745	#102-1810
Regular	#102-1681	#102-1692				#102-1749	#102-1808
Fine (50µm) ●	#102-1690	#102-1698					
Fine (40µm) ●	#102-1688						

	Curettage (con't)			Beveled cylinder	Wheel	Interproximal	
Bur No. / ISO#	878K-016	878K-018	879K-016	885-012	909-039	955-008	956-010
Diameter mm	1.6	1.8	1.6	1.2	3.9	0.8	1.0
Length mm	8.0	8.0	10.0	7.5	1.3	4.0	3.0
Coarse ●	#102-1816	#102-1827	#102-1834	#102-1845	#102-1865		
Regular	#102-1813	#102-1822	#102-1831	#102-1839	#102-1861	#102-1952	#102-1964
Fine (50µm) ●				#102-1850	#102-1869		#102-1968
Fine (40µm) ●						#102-1954	